
seca 777

Instructions for use

17-10-07-651-100d_2024-10S

English – Instructions for use	3
Français – Mode d’emploi.....	34
Español – Instrucciones de uso	66
Română – Manual de utilizare	98

CONTENTS

1 About this document	4	6 Operation	25
1.1 Representation in text	4	6.1 Switching the device on and off.....	25
1.2 Representation in diagrams	4	Switching on the device	25
2 Description of device	4	Switching off the device	25
2.1 Intended use.....	4	6.2 Switching units of measurement (variant- dependent)	25
2.2 Contraindications	4	6.3 Measuring weight	25
2.3 Clinical benefit.....	5	Weighing.....	25
2.4 Patient target group.....	5	Taring off additional weight (Tare).....	26
2.5 User qualification.....	5	Permanently displaying the weight (Hold)	26
Assembly.....	5	6.4 Measure height.....	27
Operation.....	5	Measuring heights < 1.22 m.....	27
2.6 Functional description	5	Measuring heights > 1.22 m.....	28
Measuring station.....	5	6.5 Calculating and displaying body mass index (BMI)	28
USB interface	5	7 Hygiene treatment	29
3 Safety information	5	7.1 Cleaning	29
3.1 Safety information in these instructions for use.....	5	7.2 Disinfecting.....	29
3.2 Basic safety information	6	7.3 Sterilizing.....	30
Handling device.....	6	8 Function check	30
Preventing electric shock	7	9 Troubleshooting	30
Preventing injuries and infections.....	7	10 Servicing	31
Preventing damage to device.....	8	11 Technical data	31
Handling measuring results	9	11.1 General technical data	31
Handling packaging material	9	11.2 Weight measurement	32
Handling batteries and rechargeable batteries.....	9	12 Optional accessories and spare parts	32
4 Overview	10	13 Disposal	33
4.1 Measuring station	10	13.1 Disposing of the device.....	33
4.2 Display unit.....	11	13.2 Disposing of batteries and rechargeable batteries.....	33
4.3 Display	12	14 Warranty	33
4.4 Labels	13	15 Declaration of conformity	33
5 Starting up device	15		
5.1 Scope of delivery.....	15		
5.2 Assembling device	16		
Tools required.....	16		
Fitting the measuring rod	16		
Fitting the column (white)	17		
Fitting the column (anthracite).....	18		
Fitting the display head	19		
Switching the positions of the display unit and the battery compartment.....	20		
5.3 Disinfecting the device	22		
5.4 Transporting and setting up the device.....	22		
Transporting the device.....	22		
Setting up the device.....	23		
5.5 Establishing power supply	24		
Inserting batteries.....	24		
Connecting the power supply unit (optional).....	24		

1 ABOUT THIS DOCUMENT

1.1 Representation in text

Symbol	Description
✓	Requirement for actions
▶	Action
1. 2.	Actions with specified sequence
a) b)	Steps of an action with specified sequence
⇒	Result of an action
• •	First level of a list
– –	Second level of a list

1.2 Representation in diagrams

Symbol	Description
	Indicates relevant points on the device or on device components
	Indicates directions of movement of the device or of device components
	Navigation path in menu displays

2 DESCRIPTION OF DEVICE

2.1 Intended use

The electronic measuring station supports physicians in decision-making regarding diagnosis or therapy, based on the patient's weight, height or changes of the measured parameters.

2.2 Contraindications

No contraindications are known.

2.3 Clinical benefit

The electronic measuring station supports physicians in decision-making regarding diagnosis or therapy, based on the patient's weight or height.

By that, the electronic measuring station only provides an indirect clinical benefit. The physician must perform additional examinations for an accurate diagnosis.

2.4 Patient target group

The device is intended for persons of any gender and age, with the exception of babies. There are restrictions with regard to weight and height (refer to information in the "Technical data" section).

2.5 User qualification

Assembly Devices shipped partly assembled may only be assembled by sufficiently qualified persons (such as specialist dealers, hospital technicians or seca Service technicians).

Operation The device may be operated only by medical specialists.

2.6 Functional description

Measuring station On the **seca 777** measuring station, weight is recorded by four load cells. Height is recorded with the integrated measuring rod.
Body mass index (BMI) can be calculated automatically from the measuring results.

USB interface The USB interface allows the **seca 777** measuring station to be connected to a PC by means of a suitable USB cable (USB 2.0, mini USB type B, max. cable length: 1.8 m). This makes it possible to process the weight values displayed using the customer's own software solution.

3 SAFETY INFORMATION

3.1 Safety information in these instructions for use

DANGER!

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries will occur.

WARNING!

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries may result.

CAUTION!

Used to identify a hazardous situation. If you fail to take note of this information, minor to moderate injuries may result.

NOTICE!

Used to identify possible incorrect usage of the device. If you fail to take note of this information, you may damage the device, or the measuring results may be incorrect.

NOTE

Includes additional information about use of the device.

3.2 Basic safety information

Handling device

- ▶ Please take note of the information in these instructions for use.
- ▶ Keep the instructions for use in a safe place. The instructions for use are a component of the device and must be available at all times.
- ▶ In the interest of patient safety, you and your patients are obliged to report serious events that occur in connection with this product to the manufacturer and to the authority responsible in your country.

**DANGER!****Risk of explosion**

- ▶ Do not use the device in an environment in which one of the following gases has accumulated:
 - oxygen
 - flammable anesthetics
 - other flammable substances/air mixtures

**CAUTION!****Patient hazard, damage to device**

- ▶ Additional devices connected to medical electrical devices must provide evidence of compliance with the relevant IEC or ISO standards (e.g. IEC 60950 for data-processing devices). Furthermore, all configurations must comply with the requirements of standards for medical systems (see IEC 60601-1-1 or Section 16 of edition 3.1 of IEC 60601-1 respectively). Anyone connecting additional devices to medical electrical devices is considered a system configurator and therefore responsible for ensuring that the system complies with the requirements of standards for systems. This also applies to additional devices recommended by seca. Your attention is drawn to the fact that local laws take precedence over the above-mentioned requirements of standards. In the event of any queries, please contact your local specialist dealer or Technical Service.

**CAUTION!****Patient hazard, damage to device**

- ▶ Have servicing carried out regularly as described in the relevant section of this document.
- ▶ Technical modifications may not be made to the device. The device does not contain any user-serviceable parts. Only have servicing and repairs performed by an authorized seca service partner. You can find your local service partner at www.seca.com.
- ▶ Use only seca original accessories and spare parts, otherwise seca will not grant any warranty.

**CAUTION!****Patient hazard, malfunction**

- ▶ Keep other medical electrical devices, e.g. high-frequency surgical devices, a minimum distance of approx. 1 meter away to prevent faulty measurements or wireless transmission interference.
- ▶ Keep HF devices such as cellphones at a minimum distance of approx. 1 meter to prevent faulty measurements or wireless transmission interference.
- ▶ The actual transmission output of HF devices may require minimum distances of more than 1 meter. For details, go to www.seca.com.

Preventing electric shock**WARNING!****Electric shock**

- ▶ Set up devices which can be operated with a power supply unit so that the power supply socket is within easy reach and the power supply can be disconnected quickly.
- ▶ Ensure that your local electricity supply matches the details on the power supply unit.
- ▶ Never touch the power supply unit with wet hands.
- ▶ Do not use extension cables and multiple outlets.
- ▶ Ensure that cables are not crushed or damaged by sharp edges.
- ▶ Ensure that cables do not come into contact with hot objects.
- ▶ Do not operate the device at an altitude of more than 3000 m above sea level.

Preventing injuries and infections**WARNING!****Injury from falling**

- ▶ Ensure that the device is steady and level.
- ▶ Route connector cables (if present) so that neither users nor the patient can trip over them.
- ▶ The device is not designed for supporting patients when getting up, e.g. from a wheelchair. Assist people with limited motor skills when they are getting up, e.g. from a wheelchair.
- ▶ Ensure that the patient does not step directly onto or off the edges of the weighing platform.
- ▶ Ensure that the patient steps onto and off the weighing platform slowly and safely.

**WARNING!****Danger of slipping**

- ▶ Ensure that the weighing platform is dry before the patient steps onto it.
- ▶ Ensure that the patient has dry feet before stepping onto the weighing platform.
- ▶ Ensure that the patient steps onto and off the weighing platform slowly and safely.

**WARNING!****Risk of infection**

- ▶ Before and after every measurement, wash your hands to reduce the risk of cross-contamination and nosocomial infections.
- ▶ Should the patient have contagious diseases, subject the device to a hygiene treatment as described in the relevant section of these instructions for use immediately afterwards.
- ▶ Ensure that the patient does not have any open wounds or infectious skin alterations which may come into contact with the device.
- ▶ Subject the device to a hygiene treatment at regular intervals.

Preventing damage to device**NOTICE!****Damage to device**

- ▶ Ensure that fluids never get inside the device. These can destroy the electronics.
- ▶ For devices with power supply operation: Switch off the device before you disconnect the power supply unit from the power supply socket.
- ▶ For devices with power supply operation: If the device is not to be used for an extended period, disconnect the power supply unit from the power supply socket. Only then is the device de-energized.
- ▶ For devices with battery or rechargeable battery operation: If you are not using the device for an extended period of time, remove batteries or rechargeable batteries. Only then is the device de-energized.
- ▶ Do not drop the device.
- ▶ Do not subject the device to shocks or vibrations.
- ▶ Perform a function check before each use as described in the corresponding section in this document. Do not operate the device if it is not working properly or is damaged.
- ▶ Do not place the device in direct sunlight and ensure that it is not placed in the direct proximity of a heat source. The excessive temperatures could damage the electronics.
- ▶ Avoid rapid temperature fluctuations. If the device is transported so that a temperature difference of over 20 °C occurs, the device must be left to stand for at least 2 hours before it is switched on, otherwise condensation may form; this may damage the electronics.
- ▶ Use the device only in the intended ambient conditions outlined in the section entitled "Technical data".
- ▶ Store the device only in the intended storage conditions outlined in the section entitled "Technical data".
- ▶ Use only cleaning agents and disinfectants which match the details in the section entitled "Hygiene treatment".
- ▶ For scales: Ensure that maximum capacity is not exceeded.

Handling measuring results

CAUTION! **Patient hazard**

To prevent misinterpretations, measuring results for medical purposes must only be displayed and used in SI units (kilograms/grams, meters/centimeters). Some devices have the option of displaying measuring results in different units. This is purely an additional function.

- ▶ Only use measuring results in SI units.
- ▶ The user takes sole responsibility for the use of measuring results in non-SI units.

NOTICE! **Inconsistent measuring results**

- ▶ Before you save and continue using measured values determined using this device (e.g. in seca software or in an EMR system), ensure that the measured values are plausible.
- ▶ If measured values have been sent to seca software or to an EMR system, ensure before continuing to use them that the measured values are plausible and assigned to the correct patient.

Handling packaging material

WARNING! **Risk of suffocation**

Packaging materials made of plastic film (bags) present a risk of suffocation.

- ▶ Store packaging material out of the reach of children.
- ▶ If the original packaging material is no longer available, only use plastic bags with safety holes to reduce the risk of suffocation. Use recyclable materials if possible.

NOTE
Store the original packaging material for future use (e.g. returning for servicing).

Handling batteries and rechargeable batteries

WARNING! **Injury**

Batteries and rechargeable batteries contain harmful substances which may explode and so escape if not handled properly.

- ▶ Do not try to recharge batteries.
- ▶ Do not expose batteries/rechargeable batteries to heat.
- ▶ Do not burn batteries/rechargeable batteries.
- ▶ If acid is leaking out, avoid contact with the skin, eyes, and mucous membranes. Rinse affected areas of the body with plenty of clean water and seek medical help at once.


NOTICE!

Damage to device and malfunction due to incorrect handling

- ▶ Use only the type of battery/rechargeable battery quoted in this document.
- ▶ Always replace all batteries/rechargeable batteries at the same time.
- ▶ Do not short-circuit batteries/rechargeable batteries.
- ▶ If you do not use the device for an extended period of time, remove the batteries/rechargeable batteries. This prevents acid from leaking into the device.
- ▶ If acid has penetrated the device, do not continue using it. Have the device checked by an authorized seca service partner and repaired if necessary.

4 OVERVIEW


4.1 Measuring station

Item	Device component	Function
1	Display head	<ul style="list-style-type: none">• Support element for the display unit• Housing for cable and batteries
2	Display unit	Central control and display element
3	Column	Support element for display head and measuring rod
4	Column holder	Support element for the column
5	Power supply connection socket	For connecting a seca power supply unit

Item	Device component	Function
6	Weighing platform	For weighing patients
7	Battery compartment	Takes batteries (type AA/1.5 V)
8	Measuring rod	For measuring height
9	Head slide	For sliding the measuring flap onto the lower telescopic element
10	Read-off edge 1	For heights < 1.22 m, the measuring result is read off under the head slide
11	Read-off edge 2	For heights > 1.22 m, the measuring result is read off the upper edge of the lower telescopic element
12	Lower telescopic element	For measuring height < 1.22 m
13	Upper telescopic element	For measuring height > 1.22 m
14	Measuring flap	Head stop for measuring height
15	Lock	For releasing and locking the head slide
16	Cover panel	For covering the cables
17	Recessed handle	For transporting the measuring station
18	USB interface	For transmitting the weight values displayed to the customer's software solution. If you have any questions, contact seca Service.
19	Socket for battery connector cable	Used to connect the battery connector cable
20	Socket for connector cable	Used to connect connector cable


4.2 Display unit

Item	Control	Function
1		On/off key <ul style="list-style-type: none"> Switch the device on and off
2		BMI kg/lbs key <ul style="list-style-type: none"> During the measurement procedure: <ul style="list-style-type: none"> Activate BMI mode (press key briefly) Switch between the units kilograms and pounds (press and hold key) - function depends on variant In BMI mode: <ul style="list-style-type: none"> Confirm value for height and display BMI (press key briefly) Exit BMI mode (press key briefly again)

Item	Control	Function
3		Tare ▽ key <ul style="list-style-type: none"> • During the measurement procedure: <ul style="list-style-type: none"> – Activate/deactivate Tare function • In BMI mode: <ul style="list-style-type: none"> – Reduce value for height
4		Hold △ key <ul style="list-style-type: none"> • During the measurement procedure: <ul style="list-style-type: none"> – Activate/deactivate Hold function • In BMI mode: <ul style="list-style-type: none"> – Increase value for height


4.3 Display

Item	Symbol	Meaning
A		Batteries are low
B	kg	Weight in kilograms
C	lbs	Weight in pounds
D	cm	Height in centimeters
E	BMI	Body Mass Index
F	NET	Additional weight tared off
G	HOLD	Weight is displayed permanently
H		Non-verifiable function active
I		Operation with power supply unit

4.4 Labels


Markings on the device and on the type plate	
Text/symbol	Meaning
	Name and address of manufacturer, date of manufacture
UDI	Unique Device Identification (product identification number as per Directive (EU) 2017/745)
	Article number
	Serial number
ProdID	seca product identification number
	Follow instructions for use
	Medical electrical device, Type B
d	Value in units of mass (non-verified models) States the difference between two consecutive display values
IP20	Type of protection to IEC 60529: <ul style="list-style-type: none"> • Protection against ingress of solid foreign bodies with a diameter of over 12.5 mm • Protection against access with fingers • No protection from ingress of water
	Device complies with EU directives
	Medical device in accordance with Regulation (EU) 2017/745
	Device meets the requirements of the USA and Canada. Certified and tested by a licensing laboratory (NRTL) of TÜV SÜD Product Services GmbH.
	Device complies with United Kingdom directives xxxx : Notified Body for Medical Devices of the United Kingdom
	Importer/representative in the United Kingdom: seca Ltd 40 Barn Street B5 5QB Birmingham United Kingdom
	Importer/representative in Switzerland: seca ag (schweiz) Medizinische Waagen und Messsysteme Schönmatt Str. 2 CH-4153 REINACH

Markings on the device and on the type plate	
Text/symbol	Meaning
 <p>x-y V max. xxx mA use compatible seca adapter only</p>	<p>Type plate on the power supply connection socket</p> <ul style="list-style-type: none"> • Power supply voltage required in V • Maximum current consumption in mA • : Note polarity of device connector • : Operate device with direct current • : Use compatible seca power supply units only
	Do not dispose of device in household waste
	Device can tip over. Do not push it or lean on it.

Markings on the packaging	
Symbol	Meaning
	Protect from moisture
	Arrows indicate top of product Transport and store in an upright position
	Fragile, do not throw or drop
	Permitted min. and max. temperature for transport and storage
	Permitted min. and max. humidity for transport and storage
	Permitted min. and max. air pressure for transport and storage
	Open packaging here
	Packaging material can be disposed of through recycling programs

5 STARTING UP DEVICE

5.1 Scope of delivery

Item	Component	Pcs.
a	Weighing platform with column holder	1
b	Display head with display unit and connector cable	1
c	Measuring rod	1
d	Medical variant (white column): with four preassembled screws Fitness variant (anthracite column, not shown): no preassembled screws	1
e	Spacer sleeves	2
f	Slot-head screws	2
g	Cross-head screws	2
h	Medical variant only (white column): Lock washers	4
i	Fitness variant only (anthracite column): Centering screws	4
j	T20 Allen key	1
k	Batteries (type AA/1.5 V)	6
l	Instructions for use	1

5.2 Assembling device

NOTE

- Use an assistant to perform the assembly.
- seca recommends placing the components on the floor and only setting the device in an upright position once you have fitted the measuring rod and the column.

Tools required


In addition to the scope of delivery, you require the following tools:


- Slot-head screwdriver 1x5.5 mm
- Cross-head screwdriver PH 2

Fitting the measuring rod


Proceed as follows to fit the measuring rod to the column:

1. Lie the column down on the floor.
2. Place the spacer sleeves on the bores in the column.
3. Extend the upper telescopic element far enough to expose the lower assembly bore in the lower telescopic element.
4. Position the measuring rod on the spacer sleeves as shown in the diagram.
5. Place a slot-head screw in the lower assembly bore and the spacer sleeve underneath it.
6. Tighten the slot-head screw (torque: 2 Nm).


7. Extend the upper telescopic element far enough to expose the upper assembly bore in the lower telescopic element.


8. Place a slot-head screw in the upper assembly bore and the spacer sleeve underneath it.
9. Tighten the slot-head screw (torque: 2 Nm).
10. Push the upper telescopic element back into its lowest position.

Fitting the column (white)


Proceed as follows to fit the column to the weighing platform:


1. Lie the column on the floor and position the weighing platform as shown in the diagram.
2. Hold the weighing platform steady.

3. Push the column and measuring rod into the openings in the column holder of the weighing platform.
 - ⇒ The screws on the column protrude from the underside of the weighing platform.

4. Place one lock washer in each of the bores for the screws.


5. Tighten the screws using the T20 Allen key supplied (torque: 2 Nm).

6. Set the measuring station in an upright position.


Fitting the column (anthracite)


Proceed as follows to fit the column to the weighing platform:

1. Lie the column on the floor and position the weighing platform as shown in the diagram.
2. Hold the weighing platform steady.

3. Position the column on the column holder of the weighing platform.



4. Put the centering screws in the assembly bores of the column holder and the column.
5. Tighten the screws using the T20 Allen key supplied.
6. Set the measuring station in an upright position.

Fitting the display head

Proceed as follows to fit the display head to the column:

1. Extend the upper telescopic element of the measuring rod approx. 10 cm.
2. Guide the connector cable of the display head through the center channel of the column as far as the lower part of device.
3. Put the display head on the column.
4. Push the upper telescopic element into its lowest position.


5. Place one cross-head screw in each of the two side bores of the column.
6. Tighten the cross-head screws (torque: 1 Nm).
7. Place the measuring station on the floor with the display unit facing upward.

NOTICE!

Damage due to error in installation

If the cables are bent at too sharp an angle, they may be damaged and thus cause a functional failure.

- Route all cables to prevent sharp bends and kinked connectors.



8. Connect the connector cable:
 - a) Pull the connector cable out of the column
 - b) Route the connector cable past the cable clip and plug into the socket of the weighing platform
 - c) Fix the connector cable in position using the cable clip
9. Set the measuring station in an upright position.
10. Remove the protective film from the display (if desired).


Switching the positions of the display unit and the battery compartment

At the factory, the display unit is fitted to the right-hand operating side of the display head. The battery compartment is on the left-hand operating side.


Proceed as follows to switch the positions of the display unit and the battery compartment:


1. Remove the cover panel:
 - a) Unscrew two cross-head screws
 - b) Remove the cover panel


2. Remove the battery compartment lid:
 - c) Hold the battery compartment lid steady
 - d) Remove the cross-head screw behind the battery compartment lid
 - e) Remove the battery compartment lid

3. Remove the display unit:
 - f) Undo the cross-head screw behind the display unit
 - g) Remove cable tie on the cable hook, if present
 - h) Draw cables out of the cable hook
 - i) Lift off the display unit


4. Switch the positions of the display unit and the battery holder:
 - j) Put the battery holder in its new position
 - k) Carefully put the display unit in its new position


5. Route the cables so that they do not dangle out of the display head.

NOTE


When fitting on the right-hand operating side: You can gather all the cables into the cable hook and secure them with a cable tie if necessary.

6. Mount the display unit:
 - l) Hold the display unit steady
 - m) Tighten the cross-head screw behind the display unit

7. Fit the battery compartment lid:
 - n) Insert the battery compartment lid and hold it steady
 - o) Tighten the cross-head screw behind the battery compartment lid


8. Fit the cover panel:
 - p) Insert cover panel
 - q) Place and tighten two cross-head screws

5.3 Disinfecting the device

- ▶ Disinfect the following components **before** initial commissioning as described in the section entitled “Disinfecting”:
 - Column and recessed handle
 - Display head and display unit
 - Measuring rod, head slide with measuring flap
 - Weighing platform
 - Handrail (optional accessory)

5.4 Transporting and setting up the device

Transporting the device

The device has two casters that facilitate transport over short distances.

1. If necessary, disconnect the power supply unit from the power socket and from the device.
2. Fold down the measuring flap and push the upper telescopic element of the measuring rod into its lowest position.
3. Take hold of the recessed handle of the device.

CAUTION! Injury from device falling over

- ▶ Hold onto the device firmly when transporting it.
4. Tilt the device until it can be moved freely on the casters.
 5. Move the device slowly and carefully over door sills.
 6. Transport the device to its new setup location.

Setting up the device


WARNING!

Injury from device falling over

The device is intended as a mobile medical device and is therefore not anchored permanently to a wall or the floor. A device with handrail fitted may tip over if not used in accordance with the intended purpose (e.g. as a “stand assist”).

- ▶ Ensure that the patient does not use the device as a stand assist.
- ▶ Always set up devices with a handrail at a distance of 5 to 10 cm from the wall.

The distance of at least 5 cm from the wall prevents faulty measurements. The distance of no more than 10 cm stops the device tipping over towards the wall.

CAUTION!



Faulty measurement as a result of devices being set up incorrectly

Soft floors, such as wooden boards, give under the patient's weight and falsify the measuring result. Incorrect ambient conditions or lack of alignment lead to measurement errors.

- ▶ Select a setup location with a level, stable floor in order to achieve precise measuring results.
 - ▶ Use the device only in the intended ambient conditions outlined in the section entitled “Technical data”.
 - ▶ Ensure the device is aligned before using it and every time it changes location.
- ▶ Place the device on a firm, level surface.

5.5 Establishing power supply

Inserting batteries

You need six batteries, type AA/1.5 V. Proceed as follows to establish the power supply:

1. Open the battery compartment:
 - a) Hold the battery compartment lid steady
 - b) Remove the cross-head screw behind the battery compartment lid
 - c) Remove the battery compartment lid

2. Put the new batteries in the battery holder.

NOTE

- If replacing batteries: Dispose of the old batteries as described in the section entitled “Disposal”.
- Ensure the correct polarity of the batteries (markings on the battery holder). If **bAtt** appears in the display, you have inserted one of the batteries the wrong way round or the batteries have discharged. If batteries are inserted the wrong way round, they must be removed again immediately.

3. Insert the battery holder in the battery compartment.
4. Close the battery compartment:
 - d) Insert the battery compartment lid
 - e) Tighten the cross-head screw behind the battery compartment lid

Connecting the power supply unit (optional)

The seca power supply unit is available as an optional accessory.

WARNING!


Personal injury or damage to the device as a result of incorrect power supply units

Conventional power supply units may deliver a higher voltage than is indicated on them. The measuring device may overheat, catch fire, melt or short-circuit.

- ▶ Use only original seca power supply units with a controlled 12 V output voltage.

The connection for the power supply unit is located on the side of the column holder. To establish the power supply, proceed as follows:

1. Insert the device connector of the power supply unit into the power supply connection socket of the measuring station.
2. Plug the power supply unit into a power supply socket.


6 OPERATION

6.1 Switching the device on and off

Switching on the device

1. Ensure that the device has no load.
2. Press the  key.
 - ⇒ All the elements of the display are shown briefly.
 - ⇒ **SECA** is displayed.
 - ⇒ The device is ready for operation when **0.0** is displayed.

Switching off the device

- ▶ Press the  key.

NOTE

In battery mode, the device remains on standby for 30 seconds. If there is no input or no load is put on the weighing platform during the standby time, the device switches off automatically.

6.2 Switching units of measurement (variant-dependent)

CAUTION! Patient hazard

To prevent misinterpretations, measuring results for medical purposes must only be displayed and used in SI units (kilograms/grams, meters/centimeters). Some devices have the option of displaying measuring results in different units. This is purely an additional function.

- ▶ Only use measuring results in SI units.
- ▶ The user takes sole responsibility for the use of measuring results in non-SI units.

You can switch between the metric unit kilogram and the imperial unit pounds during the measurement procedure.

- ▶ Keep the **BMI kg/lbs** key depressed until the unit of measurement display changes.

6.3 Measuring weight

Weighing

CAUTION! Injury from incorrect position of the measuring flap

Injuries may result if the measuring flap is at patient height when folded out.

- ▶ Ensure that the measuring flap is positioned well above the patient's head before the patient steps onto the weighing platform.

**CAUTION!****Faulty measurement as a result of force shunt**

If the patient is in contact with parts other than the weighing platform during the weight measurement, the device will not measure correctly.

- ▶ Ensure that the patient is in contact only with the following parts during weight measurement: Weighing platform and **seca 477** handrail (if present).

1. Ask the patient to step onto the weighing platform.
2. Ask the patient to keep still.
3. Read off the measuring result.

Taring off additional weight (Tare)

Use the **Tare** function to prevent an additional weight (e.g. a walking stick) affecting the patient's weight value.

**CAUTION!****Faulty measurement as a result of force shunt**

If an additional weight is touching the surface on which the device is standing, the weight will not be measured correctly.

- ▶ Make sure that additional weights are placed only on the weighing platform.

NOTE

The maximum weight which can be displayed is reduced by the weight of the objects which have been tared.

1. Place the additional weight on the weighing platform.
2. Press and hold the **Tare** ▾ key until the **NET** message is displayed.
3. Wait until the display unit stops flashing and **0.0** is displayed instead.

NOTE

If you remove the additional weight again, - - - - will be displayed.

4. Ask the patient to step onto the weighing platform with the additional weight (or carrying the additional weight).
 - ⇒ The patient's weight is displayed. The additional weight is deducted automatically.
5. Read off the measuring result.
6. To deactivate the **Tare** function, press and hold the **Tare** ▾ key until the **NET** message is no longer displayed.

NOTE

When you switch the device off, the **Tare** function is deactivated automatically.

Permanently displaying the weight (Hold)

If you activate the **Hold** function, the weight value is displayed permanently.

1. Press the **Hold** Δ key.
 - ⇒ The display flashes until a stable weight is measured.
 - ⇒ The weight is displayed permanently.
 - ⇒ The message Δ **HOLD** is displayed.
2. Read off the measuring result.
3. To deactivate the **Hold** function, press the **Hold** Δ key.

6.4 Measure height

You can measure height using the integrated measuring rod.


CAUTION!

Injury from incorrect position of the measuring flap


Injuries may result if the measuring flap is at patient height when folded out.

- ▶ Ensure that the measuring flap is positioned well above the patient's head before the patient steps onto the weighing platform.


Measuring heights < 1.22 m

1. Ensure that the upper telescopic element is in its lowest position.
2. Release the head slide by pressing on the lock and moving the head slide downward.


3. Fold up the measuring flap until it engages in a horizontal position.
4. Ask the patient to stand under the measuring flap and to stand up straight.
5. Push the head slide down until the measuring flap is in contact with the patient's head.


6. Read off the measuring result below the head slide (read-off edge 1) (in this case: 108.2 cm).
 - ⇒ You have the following options for continuing:
 - ▶ Make a note of the measuring result.
 - ▶ Enter the measuring result on the display unit (→ [Calculating and displaying body mass index \(BMI\)](#), page 28).
7. Ask the patient to step off the weighing platform.
8. Fold down the measuring flap.
9. Push the head slide until it engages in the lock.

Measuring heights > 1.22 m

1. Fold up the measuring flap until it engages in a horizontal position.
2. Extend the upper telescopic element until the patient can get under the measuring flap comfortably.
3. Ask the patient to stand under the measuring flap and to stand up straight.

4. Push the upper telescopic element down until the measuring flap is in contact with the patient's head.
5. Read off the measuring result on the lower telescopic element (read-off edge 2) (in this case: 145.7 cm).
⇒ You have the following options for continuing:
 - ▶ Make a note of the measuring result.
 - ▶ Enter the measuring result on the display unit (→ [Calculating and displaying body mass index \(BMI\)](#), page 28).
6. Ask the patient to step off the weighing platform.
7. Push the upper telescopic element into its lowest position.
8. Fold down the measuring flap.

6.5 Calculating and displaying body mass index (BMI)

Body Mass Index correlates height and weight.

- ✓ Weight measured
- ✓ Height measured

1. Press the **BMI kg/lbs** key briefly.
⇒ The last height entered in cm and **BMI** are displayed.
2. Use the **Tare** ▾ and **Hold** ▲ keys to enter the height in centimeters.

3. Press the **BMI kg/lbs** key briefly.
⇒ Body Mass Index is calculated and displayed automatically.
4. To deactivate the function, press the **BMI kg/lbs** key briefly.

7 HYGIENE TREATMENT

WARNING! **Electric shock**

The device is not de-energized when the on/off button is pressed and the display goes out. Use of fluids on the device may cause an electric shock.

- ▶ Before each hygiene treatment, ensure that the device is switched off.
- ▶ Disconnect the power supply connector before each hygiene treatment.
- ▶ Take the batteries out of the device before every hygiene treatment.
- ▶ Ensure that no fluids penetrate the device.

WARNING! **Risk of infection**

- ▶ Subject the device to a hygiene treatment at regular intervals as described in this section.

NOTICE! **Damage to device**

Unsuitable cleaning agents and disinfectants may damage the sensitive surfaces of the device and lead to transparent components becoming clouded.

- ▶ Do not use aggressive or abrasive cleaning agents.
- ▶ Do not use organic solvents (e.g. white spirit or petroleum spirit).

7.1 Cleaning

- ▶ If required, moisten a soft cloth with a mild soap solution and wipe the device over with it.

7.2 Disinfecting

1. Disinfect the device at regular intervals with a disinfectant suitable for sensitive surfaces and acrylic glass (e.g. 70 % ethanol).
2. Follow the instructions for use for the disinfectant.
3. Disinfect the device as follows:
 - ▶ Moisten a soft cloth with disinfectant and wipe down the device with it.
 - ▶ Comply with the intervals, see table:

Interval	Component
Before and after each measurement	<ul style="list-style-type: none"> • Weighing platform • Head slide with measuring flap • Handrail (if present)
As required	<ul style="list-style-type: none"> • Column, display unit and display head • Telescopic elements of the measuring rod • Recessed handle

7.3 Sterilizing

The device must not be sterilized.

8 FUNCTION CHECK

► Perform a function check before each use.

A complete function check includes:

- A visual inspection for mechanical damage
- A test of the alignment of the device
- A visual and function check of the display elements
- A function check of all the controls shown in the section entitled “Overview”
- A function check of the optional accessories

If you find faults or deviations during the function check, first try to remedy the fault with the aid of the “Troubleshooting” section in this document.

CAUTION! Personal injury

If you find faults or deviations during the function check which you are unable to remedy with the aid of the “Troubleshooting” section in this document, you must not use the device.

- Have the device repaired by seca Service or an authorized service partner.
- Follow the section entitled “Servicing” in this document.

9 TROUBLESHOOTING

Fault	Cause	Remedy
No weight value is displayed even though there is a load on the scale.	The device has no power supply.	<ul style="list-style-type: none"> • Switch on device. • In battery mode: Insert batteries. • In power supply operation: Connect the power supply unit.
0.0 is not displayed before weighing takes place.	A load was already on the device before it was switched on.	<ul style="list-style-type: none"> • Remove the load from the device. • Switch off the device, then switch it back on again.
----- is displayed before weighing takes place.	A load was already on the device before it was switched on.	<ul style="list-style-type: none"> • Remove the load from the device. • Switch off the device, then switch it back on again.
The weight value is flashing.	Consistent weight value not yet detected.	<ul style="list-style-type: none"> • One of the following functions is active: Hold, Autohold, Tare. Wait until the device has detected a consistent weight value. The display stops flashing. • Take the load off the device, wait until 0.0 is displayed and weigh again.
One display segment is displayed continuously or not at all.	The display is defective.	<ul style="list-style-type: none"> • Inform seca Service.
 is displayed.	Battery voltage is dropping.	<ul style="list-style-type: none"> • Replace the batteries.
bAtt is displayed.	The batteries have discharged	<ul style="list-style-type: none"> • Replace the batteries.
StOP is displayed.	Maximum capacity has been exceeded.	<ul style="list-style-type: none"> • Remove the load from the device.

Fault	Cause	Remedy
Er:[No.]:11 is displayed.	The device has too high a load or too high a load in one corner.	<ul style="list-style-type: none"> Remove the load from the device or distribute the weight more evenly. Switch off the device, then switch it back on again. If the problem persists, inform seca Service.
Er:[No.]:12 is displayed.	The device has been switched on with too high a load.	<ul style="list-style-type: none"> Remove the load from the device. Switch off the device, then switch it back on again.
Er:[No.]:16 is displayed.	The device was caused to oscillate, the zero point could not be determined.	<ul style="list-style-type: none"> Switch off the device, then switch it back on again.
Er:[No.]:32 is displayed.	An internal communication error occurred.	<ul style="list-style-type: none"> Switch off the device, then switch it back on again. If the problem persists, inform seca Service.

10 SERVICING

The product needs to be set up carefully and serviced regularly. Depending on how frequently the product is used, seca recommends servicing at intervals of 3 to 5 years.

CAUTION!

Faulty measurements as a result of poor servicing

- ▶ Have servicing and repairs carried out exclusively by seca Service or an authorized service partner.
- ▶ You can find your local service partner at www.seca.com.

11 TECHNICAL DATA

11.1 General technical data

Dimensions	
<ul style="list-style-type: none"> Depth Width Height (measuring rod not extended) 	<p>576 mm 410 mm 1356 mm</p>
Weight of measuring station	Approx. 13 kg
Ambient conditions, operation	
<ul style="list-style-type: none"> Temperature Air pressure Humidity 	<p>+10 °C to +40 °C / +50 °F to 104 °F 700 – 1060 hPa 30 % – 80 %, no condensation</p>
Ambient conditions, storage	
<ul style="list-style-type: none"> Temperature Air pressure Humidity 	<p>-10 °C to +65 °C / +14 °F to 149 °F 700 – 1060 hPa 0 % – 95 %, no condensation</p>
Ambient conditions, transport	
<ul style="list-style-type: none"> Temperature Air pressure Humidity 	<p>-10 °C to +65 °C / +14 °F to 149 °F 700 – 1060 hPa 0 % – 95 %, no condensation</p>

Power supply <ul style="list-style-type: none"> Batteries <ul style="list-style-type: none"> Supply voltage Battery type Power supply unit (optional accessory) <ul style="list-style-type: none"> Supply voltage Maximum current consumption Insulated device, protection class II (EN 60601-1) 	<p style="text-align: right;">9 V 6 x type AA/1.5 V</p> <p style="text-align: right;">12 V Typically 500 mA</p> <p style="text-align: right;"><input type="checkbox"/></p>
Mains voltage	100 V – 240 V
Power supply frequency	50 Hz – 60 Hz
Current consumption	Approx. 37 mA
Maximum runtime in battery operation (USB interface not active)	Approx. 32 hours
Metrology data, measuring rod <ul style="list-style-type: none"> Measuring range 1 Measuring range 2 Graduation Accuracy 	<p>10 cm – 122 cm (4 inch – 48 inch)</p> <p>122 cm – 230 cm (48 inch – 90.5 inch)</p> <p>1 mm (1/8 inch)</p> <p>± 5 mm</p>
Medical device in line with Directive 93/42/EEC	Class I with measuring function
EN 60601-1: Medical electrical equipment, type B	
Type of protection	IP20
Duty cycle	Continuous duty
Application parts in line with IEC 60601-1	Weighing platform, measuring flap

11.2 Weight measurement

Maximum capacity	250 kg/550 lbs
Minimum capacity	0.2 kg/0.4 lbs
Graduation	0.1 kg/0.2 lbs
Tare range	up to 250 kg (subtractive)
Accuracy <ul style="list-style-type: none"> 0 kg to 70 kg (0 lbs to 132 lbs) 70 kg to 250 kg (132 lbs to 551 lbs) 	<p>± 100 g (0.22 lbs)</p> <p>± 0.15 %</p>

12 OPTIONAL ACCESSORIES AND SPARE PARTS

Accessory/spare part	Article number
Switch-mode power supply unit: 100–240 V~ / 50–60 Hz, 12 V=/0.5 A	68 32 10 270
seca 477 handrail	477 0000 009

13 DISPOSAL

13.1 Disposing of the device

Do not dispose of the device in your household waste. The device must be properly disposed of as electronic scrap. Follow your respective national regulations. For more information, please contact seca Service or an authorized service partner.

13.2 Disposing of batteries and rechargeable batteries

Do not dispose of used batteries and rechargeable batteries in household waste, regardless of whether they contain harmful substances or not. As a consumer, you are legally obliged to dispose of batteries and rechargeable batteries via local authority collection points or trade collection points. Only dispose of batteries and rechargeable batteries once they are completely discharged.

14 WARRANTY

There is a two-year warranty period from delivery for defects attributable to poor materials or workmanship. All movable parts, e.g. batteries, cables, power supply units, rechargeable batteries etc. are exempt. Defects which come under the warranty will be repaired for the customer free of charge against proof of purchase. Additional claims cannot be considered. Costs of transport to and from seca are the responsibility of the customer if the device is located somewhere other than the customer's headquarters. In the event of transport damage, claims under warranty can only be made if the complete original packaging was used for transport and the device was secured and fastened in it according to its originally packaged condition. You should therefore keep all packaging parts.

The warranty will be voided if the device is opened by persons not expressly authorized by seca to do so.

In the event of a warranty issue, please contact your local seca office or the dealer from whom you ordered the product.

15 DECLARATION OF CONFORMITY

seca gmbh & co. kg hereby declares that the product complies with the terms of the applicable European directives and regulations. The unabridged declaration of conformity can be found at www.seca.com.

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