
# seca 787/797

## Instructions for use

17-10-07-652-100d\_2024-10S

<b>English – Instructions for use .....</b>	<b>3</b>
<b>Français – Mode d’emploi.....</b>	<b>50</b>
<b>Español – Instrucciones de uso .....</b>	<b>97</b>

# CONTENTS

<b>1 About this document</b> .....	<b>5</b>		
1.1 Representation in text .....	5		
1.2 Representation in diagrams .....	5		
<b>2 Description of device</b> .....	<b>5</b>		
2.1 Intended use .....	5		
2.2 Contraindications .....	5		
2.3 Clinical benefit .....	6		
2.4 Patient target group.....	6		
2.5 User qualification.....	6		
Assembly.....	6		
Operation.....	6		
Administration (models with interface module) .....	6		
2.6 Functional description .....	6		
Measuring station .....	6		
USB interface (seca 797 only) .....	6		
Internal interface module (seca 797 only) .....	6		
External interface module (optional accessory) .....	6		
seca connect 103 software (optional accessory) .....	6		
<b>3 Safety information</b> .....	<b>7</b>		
3.1 Safety information in these instructions for use .....	7		
3.2 Basic safety information .....	7		
Handling device.....	7		
Preventing electric shock .....	8		
Preventing injuries and infections.....	8		
Preventing damage to device.....	10		
Handling measuring results.....	10		
Handling packaging material .....	11		
Handling batteries and rechargeable batteries.....	11		
<b>4 Overview</b> .....	<b>12</b>		
4.1 Measuring station .....	12		
4.2 Display unit.....	13		
4.3 Display .....	15		
4.4 Menu.....	16		
Navigating in the menu.....	16		
Menu structure (seca 787).....	16		
Menu structure (seca 797).....	17		
4.5 Labels .....	17		
<b>5 Starting up device</b> .....	<b>20</b>		
5.1 Scope of delivery (seca 787) .....	20		
5.2 Scope of delivery (seca 797) .....	21		
5.3 Assembling device .....	22		
Tools required.....	22		
Fitting the measuring rod .....	22		
Fitting the column.....	23		
Fitting the column.....	23		
Fitting the display head .....	24		
Switching the positions of the display unit and the battery compartment.....	25		
			Fitting the scanner bracket (seca 797) .....
		5.4	Disinfecting the device .....
		5.5	Transporting and setting up the device .....
			Transporting the device.....
			Setting up the device .....
		5.6	Establishing power supply .....
			Inserting batteries.....
			Connecting the power supply unit (variant-dependent) .....
		5.7	Connecting the barcode scanner (seca 797) .....
		5.8	Setting up the WiFi connection (seca 797)...
			Setting the WiFi function .....
			Connecting device to a WiFi network (WPS).....
			Connecting device to the WiFi network (seca connect 103) .....
		<b>6</b>	<b>Operation</b> .....
		6.1	Switching the device on and off.....
			Switching on the device .....
			Switching off the device .....
		6.2	Switching units of measurement (variant-dependent) .....
		6.3	Measuring weight .....
			Weighing.....
			Taring off additional weight (Tare).....
			Permanently displaying the weight (Hold) .....
		6.4	Measure height.....
			Measuring heights < 1.22 m.....
			Entering height .....
			Measuring heights > 1.22 m.....
		6.5	Calculating and displaying body mass index (BMI) .....
		6.6	Scanning ID(s) and transmitting data via WiFi (seca 797) .....
		6.7	Other functions (menu).....
			Activating/deactivating the Autohold function (AHOLD).....
			Setting permanent display of the weight value (HOLD).....
			Setting filtering (FIL).....
			Setting the standby time (AOFF) .....
			Activating/deactivating switching of units of measurement (UnIt) (variant-dependent) .....
			Calibrating the measuring rod (HCAL).....
			Restoring factory settings (rESET) .....
			Resetting the WiFi setting (rESET) (seca 797).....
		<b>7</b>	<b>Hygiene treatment</b> .....
		7.1	Cleaning .....
		7.2	Disinfecting.....
		7.3	Sterilizing .....
		<b>8</b>	<b>Function check</b> .....

<b>9</b>	<b>Troubleshooting</b> .....	<b>43</b>
9.1	General faults .....	43
9.2	Problems with data transmission (interface module).....	44
<b>10</b>	<b>Servicing</b> .....	<b>45</b>
<b>11</b>	<b>Technical data</b> .....	<b>45</b>
11.1	General technical data.....	45
11.2	Weight measurement .....	46
<b>12</b>	<b>Optional accessories and spare parts (seca 787)</b> .....	<b>47</b>
<b>13</b>	<b>Optional accessories and spare parts (seca 797)</b> .....	<b>47</b>
<b>14</b>	<b>Compatible seca products</b> .....	<b>47</b>
<b>15</b>	<b>Disposal</b> .....	<b>48</b>
15.1	Disposing of the device.....	48
15.2	Disposing of batteries and rechargeable batteries.....	48
<b>16</b>	<b>Warranty</b> .....	<b>48</b>
<b>17</b>	<b>Declarations of conformity</b> .....	<b>48</b>
17.1	Europe .....	48
17.2	USA and Canada.....	49

# 1 ABOUT THIS DOCUMENT

## 1.1 Representation in text

Symbol	Description
✓	Requirement for actions
▶	Action
1. 2.	Actions with specified sequence
a) b)	Steps of an action with specified sequence
⇒	Result of an action
• •	First level of a list
– –	Second level of a list

## 1.2 Representation in diagrams

Symbol	Description
	Indicates relevant points on the device or on device components
	Indicates directions of movement of the device or of device components
	Navigation path in menu displays

# 2 DESCRIPTION OF DEVICE

## 2.1 Intended use

The electronic measuring station supports physicians in decision-making regarding diagnosis or therapy, based on the patient's weight, height or changes of the measured parameters.

## 2.2 Contraindications

No contraindications are known.

## 2.3 Clinical benefit

---

The electronic measuring station supports physicians in decision-making regarding diagnosis or therapy, based on the patient's weight or height.

By that, the electronic measuring station only provides an indirect clinical benefit. The physician must perform additional examinations for an accurate diagnosis.

## 2.4 Patient target group

---

The device is intended for persons of any gender and age, with the exception of babies. There are restrictions with regard to weight and height (refer to information in the "Technical data" section).

## 2.5 User qualification

---

**Assembly** Devices shipped partly assembled may only be assembled by sufficiently qualified persons (such as specialist dealers, hospital technicians or seca Service technicians).

**Operation** The device may only be operated by persons with formal training in the health-care sector or in medicine.

**Administration (models with interface module)** The device may only be set up and incorporated in a network by experienced administrators or hospital technicians.

## 2.6 Functional description

---

**Measuring station** On the **seca 787** and **seca 797** measuring stations, weight is recorded by four load cells. Height is recorded with the integrated measuring rod. Heights above 1.22 m are recorded digitally. Heights below 1.22 m are read off by the operator and recorded manually.

Body mass index (BMI) can be calculated automatically from the measuring results.

**USB interface (seca 797 only)** A barcode scanner can be connected to the USB interface. This allows users and patients to be identified by means of their barcodes and measuring results assigned in the EMR system.

**Internal interface module (seca 797 only)** With the aid of the internal interface module, the **seca 797** measuring station can transmit data to the **seca connect 103** software via a WiFi connection.

**External interface module (optional accessory)** With the aid of the **seca 452** external interface module, the **seca 787** measuring station can transmit data to the **seca connect 103** software via a WiFi or LAN connection.

**seca connect 103 software (optional accessory)** The **seca connect 103** software receives measurement data from an internal or external interface module and transmits it to an EMR system.

## 3 SAFETY INFORMATION

### 3.1 Safety information in these instructions for use

#### **DANGER!**

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries will occur.

#### **WARNING!**

Used to identify an extremely hazardous situation. If you fail to take note of this information, serious irreversible or fatal injuries may result.

#### **CAUTION!**

Used to identify a hazardous situation. If you fail to take note of this information, minor to moderate injuries may result.

#### **NOTICE!**

Used to identify possible incorrect usage of the device. If you fail to take note of this information, you may damage the device, or the measuring results may be incorrect.

#### **NOTE**

Includes additional information about use of the device.

### 3.2 Basic safety information

#### **Handling device**

- ▶ Please take note of the information in these instructions for use.
- ▶ Keep the instructions for use in a safe place. The instructions for use are a component of the device and must be available at all times.
- ▶ In the interest of patient safety, you and your patients are obliged to report serious events that occur in connection with this product to the manufacturer and to the authority responsible in your country.

#### **DANGER!**

##### **Risk of explosion**

- ▶ Do not use the device in an environment in which one of the following gases has accumulated:
  - oxygen
  - flammable anesthetics
  - other flammable substances/air mixtures

#### **CAUTION!**

##### **Patient hazard, damage to device**

- ▶ Additional devices connected to medical electrical devices must provide evidence of compliance with the relevant IEC or ISO standards (e.g. IEC 60950 for data-processing devices). Furthermore, all configurations must comply with the requirements of standards for medical systems (see IEC 60601-1-1 or Section 16 of edition 3.1 of IEC 60601-1 respectively). Anyone connecting additional devices to medical electrical devices is considered a system configurator and therefore responsible for ensuring that the system complies with the requirements of standards for systems. This also applies to additional devices recommended by seca. Your attention is drawn to the fact that local laws take precedence over the above-mentioned requirements of standards. In the event of any queries, please contact your local specialist dealer or Technical Service.

**CAUTION!****Patient hazard, damage to device**

- ▶ Have servicing carried out regularly as described in the relevant section of this document.
- ▶ Technical modifications may not be made to the device. The device does not contain any user-serviceable parts. Only have servicing and repairs performed by an authorized seca service partner. You can find your local service partner at [www.seca.com](http://www.seca.com).
- ▶ Use only seca original accessories and spare parts, otherwise seca will not grant any warranty.

**CAUTION!****Patient hazard, malfunction**

- ▶ Keep other medical electrical devices, e.g. high-frequency surgical devices, a minimum distance of approx. 1 meter away to prevent faulty measurements or wireless transmission interference.
- ▶ Keep HF devices such as cellphones at a minimum distance of approx. 1 meter to prevent faulty measurements or wireless transmission interference.
- ▶ The actual transmission output of HF devices may require minimum distances of more than 1 meter. For details, go to [www.seca.com](http://www.seca.com).

**Preventing electric shock****WARNING!****Electric shock**

- ▶ Set up devices which can be operated with a power supply unit so that the power supply socket is within easy reach and the power supply can be disconnected quickly.
- ▶ Ensure that your local electricity supply matches the details on the power supply unit.
- ▶ Never touch the power supply unit with wet hands.
- ▶ Do not use extension cables and multiple outlets.
- ▶ Ensure that cables are not crushed or damaged by sharp edges.
- ▶ Ensure that cables do not come into contact with hot objects.
- ▶ Do not operate the device at an altitude of more than 3000 m above sea level.

**Preventing injuries and infections****WARNING!****Injury from falling**

- ▶ Ensure that the device is steady and level.
- ▶ Route connector cables (if present) so that neither users nor the patient can trip over them.
- ▶ The device is not designed for supporting patients when getting up, e.g. from a wheelchair. Assist people with limited motor skills when they are getting up, e.g. from a wheelchair.
- ▶ Ensure that the patient does not step directly onto or off the edges of the weighing platform.
- ▶ Ensure that the patient steps onto and off the weighing platform slowly and safely.

**WARNING!**  
**Danger of slipping**

- ▶ Ensure that the weighing platform is dry before the patient steps onto it.
- ▶ Ensure that the patient has dry feet before stepping onto the weighing platform.
- ▶ Ensure that the patient steps onto and off the weighing platform slowly and safely.

**WARNING!**  
**Risk of infection**

- ▶ Before and after every measurement, wash your hands to reduce the risk of cross-contamination and nosocomial infections.
- ▶ Should the patient have contagious diseases, subject the device to a hygiene treatment as described in the relevant section of these instructions for use immediately afterwards.
- ▶ Ensure that the patient does not have any open wounds or infectious skin alterations which may come into contact with the device.
- ▶ Subject the device to a hygiene treatment at regular intervals.

## Preventing damage to device

### NOTICE!

#### Damage to device

- ▶ Ensure that fluids never get inside the device. These can destroy the electronics.
- ▶ For devices with power supply operation: Switch off the device before you disconnect the power supply unit from the power supply socket.
- ▶ For devices with power supply operation: If the device is not to be used for an extended period, disconnect the power supply unit from the power supply socket. Only then is the device de-energized.
- ▶ For devices with battery or rechargeable battery operation: If you are not using the device for an extended period of time, remove batteries or rechargeable batteries. Only then is the device de-energized.
- ▶ Do not drop the device.
- ▶ Do not subject the device to shocks or vibrations.
- ▶ Perform a function check before each use as described in the corresponding section in this document. Do not operate the device if it is not working properly or is damaged.
- ▶ Do not place the device in direct sunlight and ensure that it is not placed in the direct proximity of a heat source. The excessive temperatures could damage the electronics.
- ▶ Avoid rapid temperature fluctuations. If the device is transported so that a temperature difference of over 20 °C occurs, the device must be left to stand for at least 2 hours before it is switched on, otherwise condensation may form; this may damage the electronics.
- ▶ Use the device only in the intended ambient conditions outlined in the section entitled "Technical data".
- ▶ Store the device only in the intended storage conditions outlined in the section entitled "Technical data".
- ▶ Use only cleaning agents and disinfectants which match the details in the section entitled "Hygiene treatment".
- ▶ For scales: Ensure that maximum capacity is not exceeded.

## Handling measuring results

### CAUTION!

#### Patient hazard

To prevent misinterpretations, measuring results for medical purposes must only be displayed and used in SI units (kilograms/grams, meters/centimeters). Some devices have the option of displaying measuring results in different units. This is purely an additional function.

- ▶ Only use measuring results in SI units.
- ▶ The user takes sole responsibility for the use of measuring results in non-SI units.

### NOTICE!

#### Inconsistent measuring results

- ▶ Before you save and continue using measured values determined using this device (e.g. in seca software or in an EMR system), ensure that the measured values are plausible.
- ▶ If measured values have been sent to seca software or to an EMR system, ensure before continuing to use them that the measured values are plausible and assigned to the correct patient.

## Handling packaging material

### **WARNING!**

#### **Risk of suffocation**

Packaging materials made of plastic film (bags) present a risk of suffocation.

- ▶ Store packaging material out of the reach of children.
- ▶ If the original packaging material is no longer available, only use plastic bags with safety holes to reduce the risk of suffocation. Use recyclable materials if possible.

### **NOTE**

Store the original packaging material for future use (e.g. returning for servicing).

## Handling batteries and rechargeable batteries

### **WARNING!**

#### **Injury**

Batteries and rechargeable batteries contain harmful substances which may explode and so escape if not handled properly.

- ▶ Do not try to recharge batteries.
- ▶ Do not expose batteries/rechargeable batteries to heat.
- ▶ Do not burn batteries/rechargeable batteries.
- ▶ If acid is leaking out, avoid contact with the skin, eyes, and mucous membranes. Rinse affected areas of the body with plenty of clean water and seek medical help at once.


### **NOTICE!**

#### **Damage to device and malfunction due to incorrect handling**

- ▶ Use only the type of battery/rechargeable battery quoted in this document.
- ▶ Always replace all batteries/rechargeable batteries at the same time.
- ▶ Do not short-circuit batteries/rechargeable batteries.
- ▶ If you do not use the device for an extended period of time, remove the batteries/rechargeable batteries. This prevents acid from leaking into the device.
- ▶ If acid has penetrated the device, do not continue using it. Have the device checked by an authorized seca service partner and repaired if necessary.

## 4 OVERVIEW


### 4.1 Measuring station

Item	Device component	Function
1	Display head	<ul style="list-style-type: none"> <li>Support element for the display unit</li> <li>Housing for cable and batteries</li> </ul>
2	Display unit	Central control and display element
3	Column	Support element for display head and measuring rod
4	Column holder	Support element for the column
5	Power supply connection socket	For connecting a seca power supply unit
6	Weighing platform	For weighing patients
7	Battery compartment	Takes batteries (type AA/1.5 V)
8	Measuring rod	For measuring height
9	Head slide	For sliding the measuring flap onto the lower telescopic element
10	Read-off edge 1	For heights < 1.22 m, the measuring result is read off under the head slide
11	Read-off edge 2	For heights > 1.22 m, the measuring result is read off the upper edge of the lower telescopic element
12	Lower telescopic element	For measuring height < 1.22 m
13	Upper telescopic element	For measuring height > 1.22 m
14	Measuring flap	Head stop for measuring height
15	Lock	For releasing and locking the head slide

Item	Device component	Function
16	Cover panel	For covering the cables
17	Recessed handle	For transporting the measuring station
18	USB interface ( <b>seca 797</b> only)	For connecting a barcode scanner
19	Socket for electronic height measuring unit	For connecting the electronic height measuring unit
20	Socket for connector cable	Used to connect connector cable
21	Socket with blind plug	Do not use ( <b>seca 787</b> : Connection for <b>seca 452</b> ) optional accessory


## 4.2 Display unit

Item	Control	Function
1		On/off key <ul style="list-style-type: none"> <li>Switch the device on and off</li> </ul>
2		<b>Height</b> key <ul style="list-style-type: none"> <li>During the measurement procedure (measuring range 2): <ul style="list-style-type: none"> <li>Activate/deactivate permanent display of height</li> </ul> </li> <li>During the measurement procedure (measuring range 1): <ul style="list-style-type: none"> <li>Activate input of height</li> </ul> </li> <li>In the menu: <ul style="list-style-type: none"> <li>No function</li> </ul> </li> </ul>
3		<b>kg/lbs</b> key <ul style="list-style-type: none"> <li>During the measurement procedure: <ul style="list-style-type: none"> <li>Switch between kg/cm and lbs/inch (press and hold key), function depends on variant</li> </ul> </li> <li>In the menu: <ul style="list-style-type: none"> <li>No function</li> </ul> </li> </ul>

Item	Control	Function
4		<p><b>Tare</b> key</p> <ul style="list-style-type: none"> <li>• During the measurement procedure: <ul style="list-style-type: none"> <li>– Activate/deactivate <b>Tare</b> function</li> </ul> </li> <li>• In the menu: <ul style="list-style-type: none"> <li>– No function</li> </ul> </li> </ul>
5		<p><b>Menu</b> key</p> <ul style="list-style-type: none"> <li>• During the measurement procedure: <ul style="list-style-type: none"> <li>– Call up menu (press and hold key)</li> </ul> </li> <li>• In the menu: <ul style="list-style-type: none"> <li>– Go back one menu level/exit menu (press key briefly)</li> </ul> </li> </ul>
6		<p><b>△ Hold</b> key</p> <ul style="list-style-type: none"> <li>• During the measurement procedure: <ul style="list-style-type: none"> <li>– Activate/deactivate <b>Hold</b> function</li> </ul> </li> <li>• In the menu: <ul style="list-style-type: none"> <li>– Select menu item, submenu or setting</li> </ul> </li> <li>• Increase value (press and hold key to increase the value more quickly)</li> </ul>
7		<p><b>▽ BMI</b> key</p> <ul style="list-style-type: none"> <li>• If measured values are present: <ul style="list-style-type: none"> <li>– Activate <b>BMI</b> function</li> </ul> </li> <li>• In BMI mode: <ul style="list-style-type: none"> <li>– Deactivate <b>BMI</b> function</li> </ul> </li> <li>• In the menu: <ul style="list-style-type: none"> <li>– Select menu item, submenu or setting</li> <li>– Decrease value (press and hold key to decrease the value more quickly)</li> </ul> </li> </ul>
8		<p>Confirm key</p> <ul style="list-style-type: none"> <li>• If measured values are present: <ul style="list-style-type: none"> <li>– Transmit measuring results to an EMR system via <b>seca connect 103</b> software (condition: active interface module)</li> </ul> </li> <li>• In the menu: <ul style="list-style-type: none"> <li>– Confirm selected menu item</li> <li>– Save set value</li> </ul> </li> </ul>
9	Workflow LED	<p>Shows the status of data recording and transmission (condition: active interface module):</p> <ul style="list-style-type: none"> <li>• Illuminated in green: Measurement procedure is active</li> <li>• Flashing green (approx. 5 seconds): Submitting measuring results to the EMR system (depending on setting)</li> <li>• Illuminated in green (approx. 5 seconds): Measuring results successfully sent to the EMR system (depending on setting)</li> <li>• Illuminated in red (approx. 5 seconds): Error in data transmission or measurement procedure</li> </ul> <p><b>NOTE</b></p> <ul style="list-style-type: none"> <li>▶ The data recorded and transmitted are set in the <b>seca connect 103</b> software. If you have any queries, contact your administrator or hospital technician.</li> <li>▶ The <b>seca 452</b> external interface module is available as an optional accessory for the <b>seca 787</b> measuring station (→ <a href="#">Optional accessories and spare parts (seca 787)</a>, page 47).</li> </ul>

## 4.3 Display

Item	Symbol	Meaning
A		Indicates the status of the WiFi connection and the field strength of the WiFi network (condition: active interface module): <ul style="list-style-type: none"> <li>When setting up the WiFi network: <ul style="list-style-type: none"> <li>Flashing: Device connecting to a WiFi router</li> <li>Illuminated: Device and WiFi router are connected and the WiFi connection has been made</li> </ul> </li> <li>In operation: <ul style="list-style-type: none"> <li>Illuminated: The number of bars indicates the field strength of the WiFi network</li> </ul> </li> </ul>
B	<b>AP</b>	Access Point function active (service function)
C	<b>kg</b>	Weight in kilograms
D	<b>lbs</b>	Weight in pounds
E	<b>BMI</b>	Body Mass Index
F	<b>NET</b>	Additional weight tared off
G	<b>cm</b>	Height in centimeters
H	<b>ft:in</b>	Height in feet and inches
I		Height is displayed permanently
J		Operation with power supply unit
K		Batteries are low
L	<b>HOLD</b>	Weight is displayed permanently
M		Non-verifiable function active

## 4.4 Menu

### Navigating in the menu

-  1. Press and hold the **Menu** key until a menu item is displayed.  
⇒ The last menu item selected is displayed.
-  **Hold**    **BMI** 2. Press the **△ Hold** or **▽ BMI** key until the desired menu item is displayed.
-  3. Confirm the selection by pressing the  key.  
⇒ The current setting for the menu item (or a submenu item) is displayed.
-  **Hold**    **BMI** 4. To change the setting, press the **△ Hold** or **▽ BMI** key until the desired setting is displayed.

#### NOTE

To increase or reduce a numerical value more quickly, you can press and hold the **△ Hold** or **▽ BMI** key.


5. Press the  key.  
⇒ The desired setting is adopted.  
⇒ The menu closes.

#### NOTE

- To switch to the menu above without changing a setting or to exit the menu, briefly press the **Menu** key.
- If no key is pressed for a few seconds, the menu will automatically close.

### Menu structure (seca 787)


**HCAL**    Height Calibration  
**AOFF**    Auto OFF  
**AHOLd**    Auto Hold  
**FIL**        Filter

\* variant-dependent


Menu structure (seca 797)

## 4.5 Labels

Markings on the device and the type plate	
Text/symbol	Meaning
	Name and address of manufacturer, date of manufacture
<b>UDI</b>	Unique Device Identification (product identification number as per Directive (EU) 2017/745)
<b>REF</b>	Article number
<b>SN</b>	Serial number
<b>GAL</b>	Value in m/s <sup>2</sup> (model-dependent) <ul style="list-style-type: none"> <li>• Gravitational acceleration on earth</li> <li>• Depends on the intended location</li> </ul>
<b>ProdID</b>	seca product identification number

17-10-07-662-100d\_2024-10S


Markings on the device and the type plate	
Text/symbol	Meaning
	Follow instructions for use
	Medical electrical device, Type B
d	Value in units of mass (non-verified models) States the difference between two consecutive display values
IP20	Type of protection to IEC 60529: <ul style="list-style-type: none"> <li>• Protection against ingress of solid foreign bodies with a diameter of over 12.5 mm</li> <li>• Protection against access with fingers</li> <li>• No protection from ingress of water</li> </ul>
	Device complies with EU directives
	Medical device in accordance with Regulation (EU) 2017/745
	Device meets the requirements of the USA and Canada. Certified and tested by a licensing laboratory (NRTL) of TÜV SÜD Product Services GmbH.
	Symbol of the US authority Federal Communications Commission (FCC)
FCC ID	Device license number from the US Federal Communications Commission (FCC)
IC ID	Device license number from Industry Canada
	Device complies with United Kingdom directives <b>XXXX</b> : Notified Body for Medical Devices of the United Kingdom
	Importer/representative in the United Kingdom: seca Ltd 40 Barn Street B5 5QB Birmingham United Kingdom
	Importer/representative in Switzerland: seca ag (schweiz) Medizinische Waagen und Messsysteme Schönmatt Str. 2 CH-4153 REINACH
	Type plate on the power supply connection socket <ul style="list-style-type: none"> <li>• Power supply voltage required in V</li> <li>• Maximum current consumption in mA</li> <li>• : Note polarity of device connector</li> <li>• : Operate device with direct current</li> <li>• : Use compatible seca power supply units only</li> </ul>
display cable	Connection for display cable, for model <b>seca 797</b> only
	Do not dispose of device in household waste

Markings on the device and the type plate	
Text/symbol	Meaning
	Device can tip over. Do not push it or lean on it.

Markings on the packaging	
Symbol	Meaning
	Protect from moisture
	Arrows indicate top of product Transport and store in an upright position
	Fragile, do not throw or drop
	Permitted min. and max. temperature for transport and storage
	Permitted min. and max. humidity for transport and storage
	Permitted min. and max. air pressure for transport and storage
	Open packaging here
	Packaging material can be disposed of through recycling programs


## 5 STARTING UP DEVICE

### 5.1 Scope of delivery (seca 787)

Item	Component	Pcs.
a	Weighing platform with column holder	1
b	Display head with display unit and connector cable	1
c	Measuring rod	1
d	Column (with four pre-fitted screws)	1
e	Spacer sleeves	2
f	Slot-head screws	2
g	Cross-head screws	2
h	Lock washers	4
i	T20 Allen key	1
j	Batteries (type AA/1.5 V)	6
k	Instructions for use	1

## 5.2 Scope of delivery (seca 797)

Item	Component	Pcs.
a	Weighing platform with column holder	1
b	Display head with display unit, display cable and connector cable	1
c	Scanner bracket	1
d	Measuring rod	1
e	Column (with four pre-fitted screws)	1
f	Spacer sleeves	2
g	Slot-head screws	2
h	Cross-head screws	3
i	Lock washers	4
j	T20 Allen key	1
k	Power supply unit	1
l	Instructions for use	1

## 5.3 Assembling device

### NOTE

- Use an assistant to perform the assembly.
- seca recommends placing the components on the floor and only setting the device in an upright position once you have fitted the measuring rod and the column.

### Tools required


In addition to the scope of delivery, you require the following tools:


- Slot-head screwdriver 1x5.5 mm
- Cross-head screwdriver PH 2

### Fitting the measuring rod


Proceed as follows to fit the measuring rod to the column:

1. Lie the column down on the floor.
2. Place the spacer sleeves on the bores in the column.
3. Extend the upper telescopic element far enough to expose the lower assembly bore in the lower telescopic element.
4. Position the measuring rod on the spacer sleeves as shown in the diagram.
5. Place a slot-head screw in the lower assembly bore and the spacer sleeve underneath it.
6. Tighten the slot-head screw (torque: 2 Nm).


7. Extend the upper telescopic element far enough to expose the upper assembly bore in the lower telescopic element.


8. Place a slot-head screw in the upper assembly bore and the spacer sleeve underneath it.
9. Tighten the slot-head screw (torque: 2 Nm).
10. Push the upper telescopic element back into its lowest position.

### Fitting the column


Proceed as follows to fit the column to the weighing platform:


1. Lie the column on the floor and position the weighing platform as shown in the diagram.
2. Hold the weighing platform steady.

3. Push the column and measuring rod into the openings in the column holder of the weighing platform.
  - ⇒ The screws on the column protrude from the underside of the weighing platform.

4. Place one lock washer in each of the bores for the screws.


5. Tighten the screws using the T20 Allen key supplied (torque: 2 Nm).

6. Set the measuring station in an upright position.

### Fitting the display head

Proceed as follows to fit the display head to the column:


1. Extend the upper telescopic element of the measuring rod approx. 10 cm.
2. Guide the cable(s) of the display head through the center channel of the column as far as the lower part of the device.

#### NOTE

One or two cables will be present, depending on the model:

- Connector cable (**seca 787**)
- Connector cable and display cable (**seca 797**)

3. Put the display head on the column.
4. Push the upper telescopic element into its lowest position.

5. Place one cross-head screw in each of the two side bores of the column.
6. Tighten the cross-head screws (torque: 1 Nm).
7. Place the measuring station on the floor with the display unit facing upward.


#### NOTICE!

##### Damage due to error in installation

If the cables are bent at too sharp an angle, they may be damaged and thus cause a functional failure.

- Route all cables to prevent sharp bends and kinked connectors.



8. Connect the connector cable:
  - a) Pull the connector cable out of the column
  - b) Route the connector cable past the cable clip and plug into the socket of the weighing platform
  - c) Fix the connector cable in position using the cable clip
9. **seca 797** only: Connect the display cable:
  - d) Draw the display cable out of the column
  - e) Plug the display cable into the labeled socket (**display cable**)
  - f) Fix the display cable in position with the cable clip
10. Set the measuring station in an upright position.
11. Remove the protective film from the display (if desired).


### Switching the positions of the display unit and the battery compartment

At the factory, the display unit is fitted to the right-hand operating side of the display head. The battery compartment is on the left-hand operating side.


Proceed as follows to switch the positions of the display unit and the battery compartment:


1. Remove the cover panel:
  - a) Unscrew two cross-head screws
  - b) Remove the cover panel


2. Remove the battery compartment lid:
  - c) Hold the battery compartment lid steady
  - d) Remove the cross-head screw behind the battery compartment lid
  - e) Remove the battery compartment lid

3. Remove the display unit:
  - f) Undo the cross-head screw behind the display unit
  - g) Remove both cable tie on cable hook and the strain relief, if present
  - h) Draw cables out of the cable hook
  - i) Lift off the display unit

4. Switch the positions of the display unit and the battery holder:
  - j) Put the battery holder in its new position
  - k) Carefully put the display unit in its new position



5. Route the cables so that they do not dangle out of the display head.


**NOTE**

- When fitting on the right-hand operating side: You can gather all the cables into the cable hook and secure them with a cable tie if necessary.
- **seca 797** only: If a barcode scanner is installed, route the barcode scanner cable over the strain relief and secure with a cable tie if required (→ [Connecting the barcode scanner \(seca 797\)](#), page 31).

6. Mount the display unit:
  - l) Hold the display unit steady
  - m) Tighten the cross-head screw behind the display unit


7. Fit the battery compartment lid:
  - n) Insert the battery compartment lid and hold it steady
  - o) Tighten the cross-head screw behind the battery compartment lid

8. Fit the cover panel:
  - p) Insert cover panel
  - q) Place and tighten two cross-head screws

### Fitting the scanner bracket (seca 797)


The scanner bracket can be fitted to either side of the display head. It can be fitted offset by 90° to point forwards, backwards, or sideways.

#### NOTE

If a **seca 477** handrail (optional accessory) is to be fitted, the scanner bracket must be fitted pointing sideways.

1. Place the scanner bracket in the desired position.
2. Fit the scanner bracket to the display head using a cross-head screw.


## 5.4 Disinfecting the device

- ▶ Disinfect the following components **before** initial commissioning as described in the section entitled "Disinfecting":
  - Column and recessed handle
  - Display head and display unit
  - Measuring rod, head slide with measuring flap
  - Weighing platform
  - Handrail (optional accessory)

## 5.5 Transporting and setting up the device

### Transporting the device

The device has two casters that facilitate transport over short distances.

1. If necessary, disconnect the power supply unit from the power socket and from the device.
2. Fold down the measuring flap and push the upper telescopic element of the measuring rod into its lowest position.
3. Take hold of the recessed handle of the device.

#### **CAUTION!**

##### **Injury from device falling over**

- ▶ Hold onto the device firmly when transporting it.

4. Tilt the device until it can be moved freely on the casters.
5. Move the device slowly and carefully over door sills.
6. Transport the device to its new setup location.

### Setting up the device

#### **CAUTION!**

##### **Faulty measurement as a result of devices being set up incorrectly**

Soft floors, such as wooden boards, give under the patient's weight and falsify the measuring result. Incorrect ambient conditions or lack of alignment lead to measurement errors.


- ▶ Select a setup location with a level, stable floor in order to achieve precise measuring results.
- ▶ Use the device only in the intended ambient conditions outlined in the section entitled "Technical data".
- ▶ Ensure the device is aligned before using it and every time it changes location.

**WARNING!****Injury from device falling over**

The device is intended as a mobile medical device and is therefore not anchored permanently to a wall or the floor. A device with handrail fitted may tip over if not used in accordance with the intended purpose (e.g. as a “stand assist”).

- ▶ Ensure that the patient does not use the device as a stand assist.
- ▶ Always set up devices with a handrail at a distance of 5 to 10 cm from the wall.



The distance of at least 5 cm from the wall prevents faulty measurements. The distance of no more than 10 cm stops the device tipping over towards the wall.

- ▶ Place the device on a firm, level surface.

## 5.6 Establishing power supply

### Inserting batteries

You need six batteries, type AA/1.5 V. Proceed as follows to establish the power supply:

1. Open the battery compartment:
  - a) Hold the battery compartment lid steady
  - b) Remove the cross-head screw behind the battery compartment lid
  - c) Remove the battery compartment lid

2. Put the new batteries in the battery holder.

#### NOTE

- If replacing batteries: Dispose of the old batteries as described in the section entitled “Disposal”.
- Ensure the correct polarity of the batteries (markings on the battery holder). If **bAtt** appears in the display, you have inserted one of the batteries the wrong way round or the batteries have discharged. If batteries are inserted the wrong way round, they must be removed again immediately.

3. Insert the battery holder in the battery compartment.
4. Close the battery compartment:
  - d) Insert the battery compartment lid
  - e) Tighten the cross-head screw behind the battery compartment lid

### Connecting the power supply unit (variant-dependent)

The seca power supply unit is included in the scope of delivery or is available as an optional accessory depending on the variant.

#### WARNING!


**Personal injury or damage to the device as a result of incorrect power supply units**

Conventional power supply units may deliver a higher voltage than is indicated on them. The measuring device may overheat, catch fire, melt or short-circuit.

- ▶ Use only original seca power supply units with a controlled 12 V output voltage.

The connection for the power supply unit is located on the side of the column holder. To establish the power supply, proceed as follows:

1. Insert the device connector of the power supply unit into the power supply connection socket of the measuring station.
2. Plug the power supply unit into a power supply socket.

## 5.7 Connecting the barcode scanner (seca 797)

A barcode scanner can be connected to the USB interface of the device. The barcode scanner can be used to record the ID of the patient and/or the ID of the user.


### NOTE

- If you want to use the WiFi function or a barcode scanner (variant-dependent), the power consumption of the device will rise significantly. In these cases, operate the device with the power supply unit provided.
- Observe the maximum permitted current consumption of the barcode scanner (→ [General technical data, page 45](#)).
- Use only barcode scanners recommended by seca (→ [Optional accessories and spare parts \(seca 797\), page 47](#)).

To connect a barcode scanner to the display unit, proceed as follows:


1. Remove the cover panel:

- a) Unscrew two cross-head screws
- b) Remove the cover panel

2. Connect the barcode scanner:

- c) Plug the USB stick of the barcode scanner into the USB interface of the display unit
- d) Route the barcode scanner cable over the strain relief

### NOTE


Some barcode scanner cables are very thin and are thus not adequately held in place.

- ▶ Additionally fix the barcode scanner cable in the bores on the strain relief using a cable tie.

3. Hang the barcode scanner in the scanner bracket.

4. Fit the cover panel:

- e) Insert cover panel
- f) Place and tighten two cross-head screws

## 5.8 Setting up the WiFi connection (seca 797)

### Setting the WiFi function

1. Select the **rF** menu item (→ [Navigating in the menu, page 16](#)).
2. Confirm your selection.
3. Select the **SYS** menu item.
4. Confirm your selection.  
⇒ The current setting is displayed.
5. Select the desired setting:

Option	Effect
<b>On</b>	WiFi function is activated with power supply operation, is deactivated with battery mode
<b>OFF</b>	WiFi function deactivated
<b>PErM</b>	WiFi function is activated with power supply operation and battery mode (leads to increased power consumption with battery mode)

6. Confirm your selection.  
⇒ You exit the menu automatically.

### Connecting device to a WiFi network (WPS)

If you have access to the WiFi router, you can connect your device to the WiFi network via WPS.

1. Select the **rF** menu item (→ [Navigating in the menu, page 16](#)).
2. Activate the WPS function of your WiFi router as described in the instructions for use of your WiFi router.
3. Select the **Lrn** menu item.  
⇒ The  symbol flashes.  
⇒ The device connects to the WiFi router.  
⇒ As soon as the device is connected to the WiFi network, the  symbol is permanently on.

#### NOTE

Further settings need to be made to enable measurement data to be transmitted to an EMR system via the **seca connect 103** software.

- Observe the **seca 103/452** system instructions for use.

### Connecting device to the WiFi network (seca connect 103)

If you do not have access to the WPS function of the WiFi router or want to integrate several devices, you can connect your device to the WiFi network via the **seca connect 103** software.

- ✓ Barcode scanner connected (→ [Connecting the barcode scanner \(seca 797\), page 31](#))
1. Enter the data for the WiFi network in the **seca connect 103** software as described in the **seca 103/452** system instructions for use.
  2. Scan the QR code generated in the software.  
⇒ The  symbol flashes.  
⇒ The name of the network and the password for the WiFi network are saved in the device.  
⇒ The device connects to the WiFi router.  
⇒ As soon as the device is connected to the WiFi network, the  symbol is permanently on.

**NOTE**

Further settings need to be made to enable measurement data to be transmitted to an EMR system via the **seca connect 103** software.


- ▶ Observe the **seca 103/452** system instructions for use.

## 6 OPERATION

### 6.1 Switching the device on and off

#### Switching on the device

1. Ensure that the device has no load.
2. Press the  key.
  - ⇒ All the elements of the display are shown briefly.
  - ⇒ **SECA** is displayed.
  - ⇒ The device is ready for operation when **0.0** is displayed.

#### Switching off the device

- ▶ Press the  key.

**NOTE**

In battery mode, the device remains on standby for 30 seconds. If there is no input or no load is put on the weighing platform during the standby time, the device switches off automatically.

### 6.2 Switching units of measurement (variant-dependent)

**CAUTION!****Patient hazard**

To prevent misinterpretations, measuring results for medical purposes must only be displayed and used in SI units (kilograms/grams, meters/centimeters). Some devices have the option of displaying measuring results in different units. This is purely an additional function.

- ▶ Only use measuring results in SI units.
- ▶ The user takes sole responsibility for the use of measuring results in non-SI units.

You can switch between the metric units kilograms and centimeters and the imperial units pounds and feet/inches during the measurement procedure. The units of weight and height change at the same time.

- ▶ Press and hold the **kg/lbs** key until display of the unit of measurement changes.

## 6.3 Measuring weight

### Weighing

#### CAUTION!

##### Injury from incorrect position of the measuring flap

Injuries may result if the measuring flap is at patient height when folded out.

- ▶ Ensure that the measuring flap is positioned well above the patient's head before the patient steps onto the weighing platform.

#### CAUTION!

##### Faulty measurement as a result of force shunt

If the patient is in contact with parts other than the weighing platform during the weight measurement, the device will not measure correctly.

- ▶ Ensure that the patient is in contact only with the following parts during weight measurement: Weighing platform and **seca 477** handrail (if present).

1. Ask the patient to step onto the weighing platform.
2. Ask the patient to keep still.
3. Read off the measuring result.

### Taring off additional weight (Tare)

Use the **Tare** function to prevent an additional weight (e.g. a walking stick) affecting the patient's weight value.

#### CAUTION!

##### Faulty measurement as a result of force shunt

If an additional weight is touching the surface on which the device is standing, the weight will not be measured correctly.

- ▶ Make sure that additional weights are placed only on the weighing platform.

#### NOTE

The maximum weight which can be displayed is reduced by the weight of the objects which have been tared.

1. Place the additional weight on the weighing platform.
2. Press the **Tare** key.  
⇒ The message **NET** is displayed.
3. Wait until the display stops flashing and **0.0** is displayed instead.

#### NOTE

If you remove the additional weight again, - - - - will be displayed.

4. Ask the patient to step onto the weighing platform with the additional weight (or carrying the additional weight).  
⇒ The patient's weight is displayed. The additional weight is deducted automatically.
5. Read off the measuring result.
6. To deactivate the **Tare** function, press the **Tare** key.

#### NOTE

When you switch the device off, the **Tare** function is deactivated automatically.

### Permanently displaying the weight (Hold)

If you activate the **Hold** function, the weight value is displayed permanently.

1. Press the  $\Delta$  **Hold** key.
  - ⇒ The display flashes until a stable weight is measured.
  - ⇒ The weight is displayed permanently.
  - ⇒ The message  $\Delta$  **HOLD** is displayed.
2. Read off the measuring result.
3. To deactivate the **Hold** function, press the  $\Delta$  **Hold** key.

#### NOTE

You can select which action automatically clears the displayed weight value from the display (→ [Setting permanent display of the weight value \(HOLD\)](#), page 39).

## 6.4 Measure height

You can measure height using the integrated measuring rod.



#### CAUTION!

##### Injury from incorrect position of the measuring flap


Injuries may result if the measuring flap is at patient height when folded out.

- ▶ Ensure that the measuring flap is positioned well above the patient's head before the patient steps onto the weighing platform.

### Measuring heights < 1.22 m

1. Ensure that the upper telescopic element is in its lowest position.
2. Release the head slide by pressing on the lock and moving the head slide downward.
3. Fold up the measuring flap until it engages in a horizontal position.
4. Ask the patient to stand under the measuring flap and to stand up straight.
5. Push the head slide down until the measuring flap is in contact with the patient's head.

### Entering height

6. Read off the measuring result below the head slide (read-off edge 1).  
⇒ You have the following options for continuing:
  - ▶ Make a note of the measuring result.
  - ▶ Enter the measuring result on the display unit (→ [Entering height, page 36](#)).
7. Ask the patient to step off the weighing platform.
8. Fold down the measuring flap.
9. Push up the head slide until it engages in the lock.

If you are measuring in the range < 1.22 m, the height must be entered on the display unit in order to use the following functions:

- Calculating BMI (→ [Calculating and displaying body mass index \(BMI\), page 37](#))
- Transmitting measured values to the **seca connect 103** software via WIFI (→ [Scanning ID\(s\) and transmitting data via WIFI \(seca 797\), page 38](#))

### Measuring heights > 1.22 m


1. Determine the height (→ [Measure height, page 35](#)).
2. Press the **Height** key.  
⇒ The message **SEt** and the last height entered are displayed.
3. Use the **△ Hold** and **▽ BMI** key to enter the height.
4. Press the  key.  
⇒ The height is permanently displayed.



### NOTE

During height measurement, the value in the display changes as follows:

- The value in the display changes as long as you are moving the upper telescopic element.
- The value flashes as soon as you stop moving the telescopic element.
- After approx. 20 s, the value stops flashing and - - - - is displayed. Press the **Height** key to permanently display the value.

1. Fold up the measuring flap until it engages in a horizontal position.
2. Extend the upper telescopic element until the patient can get under the measuring flap comfortably.
3. Ask the patient to stand under the measuring flap and to stand up straight.

4. Push the upper telescopic element down until the measuring flap is in contact with the patient's head.  
⇒ The value flashes in the display.

5. Press the **Height** key.  
⇒ The height and the  symbol are permanently displayed.
6. Read off the measuring result.

**NOTE**

You can also read off the measuring result on the lower telescopic element (read-off edge 2).

7. To exit permanent display of the height, press the **Height** key.  
⇒ The value flashes in the display.
8. Ask the patient to step off the weighing platform.
9. Push the upper telescopic element into its lowest position.
10. Fold down the measuring flap.

## 6.5 Calculating and displaying body mass index (BMI)

Body Mass Index correlates height and weight.

1. Determine the patient's weight (→ [Measuring weight, page 34](#)).
2. Press the **△ Hold** key.  
⇒ The weight is displayed permanently in the display.

**NOTE**

When the **Autohold** function is active, the weight is automatically displayed permanently.

3. Determine the patient's height (→ [Measure height, page 35](#)).  
⇒ The height is displayed permanently in the display.


4. Press the **▽ BMI** key.  
⇒ Body Mass Index is calculated and displayed automatically.
5. To deactivate the **BMI** function, press the **▽ BMI** key.

**NOTE**

After approx. 10 s, the display automatically displays weight and height again.

## 6.6 Scanning ID(s) and transmitting data via WiFi (seca 797)

The internal interface module of the device wirelessly connects the measuring device to an EMR system via the **seca connect 103** software.

A barcode scanner can be used to record the ID of the patient and/or the ID of the user. A WiFi connection transmits the measuring results to the electronic patient file via the **seca connect 103** software.

The following conditions need to be met in order to use these functions:

- WiFi function of the device is active
- Device is connected to an EMR system via the **seca connect 103** software
- Barcode scanner is connected to the device

The measurement procedure consists of the following steps:

- Record ID(s) using the barcode scanner
- Record weight and/or height
- Transmit measuring results to the **seca connect 103** software

### NOTE

The height must be permanently displayed (**Height** key) for it to be transmitted to the **seca connect 103** software via WiFi. The measurement procedure can be adapted on an individual basis. The options that are active depend on the settings made in the **seca connect 103** software.

1. Observe the system instructions for use of **seca 103/452** (section “7.5 Administering seca measuring devices” and “10. Operating connected seca measuring devices”).
2. Contact your administrator or hospital technician.

## 6.7 Other functions (menu)

Other functions and settings are available in the menu. For a summary of the menu structure, see the section entitled “Overview”.

### NOTE

**seca 797** only: Explanation of the functions under the **rF** menu item: → [Setting up the WiFi connection \(seca 797\)](#), page 32.

### Activating/deactivating the Autohold function (AHOLD)

If the **Autohold** function is active, the weight value is automatically displayed permanently. It is then no longer necessary to press the **Δ Hold** key.

### NOTE

You can select which action automatically clears the displayed weight value from the display (→ [Setting permanent display of the weight value \(HOLD\)](#), page 39).

1. Select the **AHOLD** menu item.
2. Confirm your selection.  
⇒ The current setting is displayed.

3. Select the desired setting:
  - **On**
  - **OFF**
4. Confirm your selection.
  - ⇒ The setting is adopted.
  - ⇒ The menu closes.

### Setting permanent display of the weight value (HOLD)

If the **Autohold** function is activated, or you have pressed the **△ Hold** key, the weight value is displayed permanently. You can select when the weight value display is automatically cleared again.

1. From the menu, select the **HOLD** menu item.
2. Confirm your selection.
  - ⇒ The current setting is displayed.
3. Select the desired setting:
  - **HOLD1**: A new load on the measuring station clears the display value
  - **HOLD2**: Removing the load from the measuring station clears the display value
4. Confirm your selection.
  - ⇒ The setting is adopted.
  - ⇒ The menu closes.

### Setting filtering (FIL)

Filtering reduces interference (due to patient movement, for example) when measuring weight. You can modify the strength of filtering.

1. From the menu, select the **FIL** menu item.
2. Confirm your selection.
  - ⇒ The current setting is displayed.
3. Select a filtering level:
  - **0**: low filtering
  - **1**: moderate filtering
  - **2**: high filtering
4. Confirm your selection.
  - ⇒ The setting is adopted.
  - ⇒ The menu closes.

### Setting the standby time (AOFF)

In battery mode the device remains on standby for 30 seconds (factory setting). If there is no input or no load is put on the weighing platform during the standby time, the device switches off automatically. Measuring results will be discarded.

You can adjust standby time as follows:

Function	Setting
Factory setting	30 seconds
Increment	30 seconds
Minimum standby time	30 seconds
Maximum standby time	600 seconds (10 minutes)

1. From the menu, select the **AOFF** menu item.
2. Confirm your selection.
  - ⇒ The time currently set is displayed (in this case: 30 seconds).
3. Select the desired time:
  - To increase value: **△ Hold** key
  - To reduce value: **▽ BMI** key
4. Confirm your selection.
  - ⇒ The setting is adopted.
  - ⇒ The menu closes.

### Activating/deactivating switching of units of measurement (Unlt) (variant-dependent)

With some device variants, you can switch the units of measurement using the **kg/lbs** key. You can switch the function of the key on and off.


1. From the menu, select the **Unlt** menu item.
2. Confirm your selection.  
⇒ The current setting is displayed.
3. Select the desired setting:
  - **On**
  - **OFF**
4. Confirm your selection.  
⇒ The setting is adopted.  
⇒ The menu closes.

### Calibrating the measuring rod (HCAL)

Calibrate the measuring rod if the value on the display deviates from the value on read-off edge 2.

1. From the menu, select the **HCAL** menu item.
2. Confirm your selection.  
⇒ **SEt** and **150.0 cm** are displayed.
3. Set the measuring rod to 150 cm.
4. Confirm the setting.  
⇒ The measuring rod is calibrated.  
⇒ The menu closes.

### Restoring factory settings (rESEt)

If you restore the factory settings (variant-dependent), the settings for the following functions will be reset:

- Autohold (**AHOLD**)
- Permanent display of weight value (**HOLD**)
- Filtering (**FIL**)
- Standby time (**AOFF**)
- Switching units of measurement (**Unit**)

Units of height and weight are reset to metric units.

#### NOTE

**seca 797** only: If you restore the factory settings, the WiFi settings will be retained.

1. From the menu, select the **rESEt** menu item.
2. Confirm your selection.  
⇒ The menu closes.
3. Switch off the device.  
⇒ Factory settings are restored and are available when the device is switched back on.

### Resetting the WiFi setting (rESEt) (seca 797)

If you reset the WiFi settings, the following information will be cleared:

- Network name (SSID)
- Network code

#### NOTE

If you reset the WiFi settings, the factory settings for the device will also be automatically restored. The factory settings will be available once the device has been restarted.

1. From the menu, select the **rESEt** menu item.
2. Press and hold the key for approx. 6 s until **0.0** is displayed.  
⇒ The WiFi settings are reset.

## 7 HYGIENE TREATMENT

### **WARNING!** **Electric shock**

The device is not de-energized when the on/off button is pressed and the display goes out. Use of fluids on the device may cause an electric shock.

- ▶ Before each hygiene treatment, ensure that the device is switched off.
- ▶ Disconnect the power supply connector before each hygiene treatment.
- ▶ Take the batteries out of the device before every hygiene treatment.
- ▶ Ensure that no fluids penetrate the device.

### **WARNING!** **Risk of infection**

- ▶ Subject the device to a hygiene treatment at regular intervals as described in this section.

### **NOTICE!**

#### **Damage to device**

Unsuitable cleaning agents and disinfectants may damage the sensitive surfaces of the device and lead to transparent components becoming clouded.

- ▶ Do not use aggressive or abrasive cleaning agents.
- ▶ Do not use organic solvents (e.g. white spirit or petroleum spirit).

### 7.1 Cleaning

- ▶ If required, moisten a soft cloth with a mild soap solution and wipe the device over with it.

### 7.2 Disinfecting

1. Disinfect the device at regular intervals with a disinfectant suitable for sensitive surfaces and acrylic glass (e.g. 70 % ethanol).
2. Follow the instructions for use for the disinfectant.
3. Disinfect the device as follows:
  - ▶ Moisten a soft cloth with disinfectant and wipe down the device with it.
  - ▶ Comply with the intervals, see table:

Interval	Component
Before and after each measurement	<ul style="list-style-type: none"><li>• Weighing platform</li><li>• Head slide with measuring flap</li><li>• Handrail (if present)</li></ul>
As required	<ul style="list-style-type: none"><li>• Column, display unit and display head</li><li>• Telescopic elements of the measuring rod</li><li>• Recessed handle</li></ul>

## 7.3 Sterilizing

The device must not be sterilized.

## 8 FUNCTION CHECK

► Perform a function check before each use.

A complete function check includes:

- A visual inspection for mechanical damage
- A test of the alignment of the device
- A visual and function check of the display elements
- A function check of all the controls shown in the section entitled “Overview”
- A function check of the optional accessories

If you find faults or deviations during the function check, first try to remedy the fault with the aid of the “Troubleshooting” section in this document.

### CAUTION!

#### Personal injury

If you find faults or deviations during the function check which you are unable to remedy with the aid of the “Troubleshooting” section in this document, you must not use the device.

- Have the device repaired by seca Service or an authorized service partner.
- Follow the section entitled “Servicing” in this document.

## 9 TROUBLESHOOTING

### 9.1 General faults

Fault	Cause	Remedy
No weight value is displayed even though there is a load on the scale.	The device has no power supply.	<ul style="list-style-type: none"> <li>• Switch on device.</li> <li>• In battery mode: Insert batteries.</li> <li>• In power supply operation: Connect the power supply unit.</li> </ul>
<b>0.0</b> is not displayed before weighing takes place.	A load was already on the device before it was switched on.	<ul style="list-style-type: none"> <li>• Remove the load from the device.</li> <li>• Switch off the device, then switch it back on again.</li> </ul>
----- is displayed before weighing takes place.	A load was already on the device before it was switched on.	<ul style="list-style-type: none"> <li>• Remove the load from the device.</li> <li>• Switch off the device, then switch it back on again.</li> </ul>
The weight value is flashing.	Consistent weight value not yet detected.	<ul style="list-style-type: none"> <li>• One of the following functions is active: <b>Hold</b>, <b>Autohold</b>, <b>Tare</b>. Wait until the device has detected a consistent weight value. The display stops flashing.</li> <li>• Take the load off the device, wait until <b>0.0</b> is displayed and weigh again.</li> </ul>
One display segment is displayed continuously or not at all.	The display is defective.	<ul style="list-style-type: none"> <li>• Inform seca Service.</li> </ul>

Fault	Cause	Remedy
 is displayed.	Battery voltage is dropping.	<ul style="list-style-type: none"> <li>Replace the batteries.</li> </ul>
<b>bAtt</b> is displayed.	The batteries have discharged.	<ul style="list-style-type: none"> <li>Replace the batteries.</li> </ul>
<b>StOP</b> is displayed.	Maximum capacity has been exceeded.	<ul style="list-style-type: none"> <li>Remove the load from the device.</li> </ul>
<b>Er:[No.]:11</b> is displayed.	The device has too high a load or too high a load in one corner.	<ul style="list-style-type: none"> <li>Remove the load from the device or distribute the weight more evenly.</li> <li>Switch off the device, then switch it back on again.</li> <li>If the problem persists, inform seca Service.</li> </ul>
<b>Er:[No.]:12</b> is displayed.	The device has been switched on with too high a load.	<ul style="list-style-type: none"> <li>Remove the load from the device.</li> <li>Switch off the device, then switch it back on again.</li> </ul>
<b>Er:[No.]:16</b> is displayed.	The device was caused to oscillate, the zero point could not be determined.	<ul style="list-style-type: none"> <li>Switch off the device, then switch it back on again.</li> </ul>
<b>Er:[No.]:32</b> is displayed.	An internal communication error occurred.	<ul style="list-style-type: none"> <li>Switch off the device, then switch it back on again.</li> <li>Ensure that the display cable is connected correctly (<b>seca 797</b>).</li> <li>If the problem persists, inform seca Service.</li> </ul>

## 9.2 Problems with data transmission (interface module)

### NOTE

Observe the system instructions for use of **seca 103/452** and the instructions for use of the barcode scanner.

Fault	Cause	Remedy
Workflow LED does not come on	The WiFi function is deactivated	<ul style="list-style-type: none"> <li>Activate WiFi function</li> </ul>
	Workflow LED is defective	<ul style="list-style-type: none"> <li>Contact your administrator or hospital technician</li> <li>Inform seca Service</li> </ul>
The Workflow LED illuminates in red, the measuring result is displayed	Errors occurred during data transmission	<ul style="list-style-type: none"> <li>Ask the patient to step off the device and wait until the Workflow LED goes out</li> <li>Check WiFi connection</li> <li>Deactivate <b>BMIF</b> function</li> <li>Repeat measurement procedure</li> <li>If the problem persists, contact your administrator or hospital technician</li> </ul>
The Workflow LED illuminates in red, an error code ( <b>Er:8:[No.]</b> ) is displayed	Errors occurred during the measurement procedure or data transmission	<ul style="list-style-type: none"> <li>Contact your administrator or hospital technician</li> </ul>
<b>Er:8:94</b> is displayed	An error occurred during the measurement procedure: IDs were scanned at the wrong time	<ul style="list-style-type: none"> <li>Cancel measurement procedure</li> <li>Repeat measurement procedure and observe proper sequence</li> <li>If the problem persists, contact your administrator or hospital technician</li> </ul>
<b>Er:8:95</b> is displayed	Hardware fault	<ul style="list-style-type: none"> <li>Contact your administrator or hospital technician</li> <li>Inform seca Service</li> </ul>
<b>Er:[No.]:71</b> is displayed after the  key is pressed	The WiFi function is deactivated	Activate WiFi function

Fault	Cause	Remedy
The <b>rF</b> menu item is not visible	The internal interface module of the device is defective	<ul style="list-style-type: none"> <li>• Inform seca Service</li> </ul>
Only the <b>SYS</b> menu item is visible in the <b>rF</b> menu	The WiFi function is deactivated	Activate WiFi function
<b>noChG</b> is displayed if the WiFi function activated ( <b>rF\SYS\On</b> )	The device is in battery mode	<ul style="list-style-type: none"> <li>• Use power supply unit (recommended)</li> <li>• In the <b>rF\SYS</b> menu, select the option <b>PErM</b></li> </ul>
<b>Id:Er</b> is displayed	The ID of the user or the patient has not been found in the EMR system	<ul style="list-style-type: none"> <li>• Rescan ID</li> <li>• If the problem persists, contact your administrator or hospital technician</li> </ul>

## 10 SERVICING

The product needs to be set up carefully and serviced regularly. Depending on how frequently the product is used, seca recommends servicing at intervals of 3 to 5 years.

### CAUTION!

Faulty measurements as a result of poor servicing

- ▶ Have servicing and repairs carried out exclusively by seca Service or an authorized service partner.
- ▶ You can find your local service partner at [www.seca.com](http://www.seca.com).

## 11 TECHNICAL DATA

### 11.1 General technical data

Dimensions <ul style="list-style-type: none"> <li>• Depth</li> <li>• Width</li> <li>• Height (measuring rod not extended)</li> </ul>	<p>576 mm</p> <p>410 mm</p> <p>1356 mm</p>
Weight of measuring station <ul style="list-style-type: none"> <li>• <b>seca 787</b></li> <li>• <b>seca 797</b></li> </ul>	<p>Approx. 12.9 kg</p> <p>Approx. 13.1 kg</p>
Ambient conditions, operation <ul style="list-style-type: none"> <li>• Temperature</li> <li>• Air pressure</li> <li>• Humidity</li> </ul>	<p>+10 °C to +40 °C / +50 °F to 104 °F</p> <p>700 – 1060 hPa</p> <p>30 % – 80 %, no condensation</p>
Ambient conditions, storage <ul style="list-style-type: none"> <li>• Temperature</li> <li>• Air pressure</li> <li>• Humidity</li> </ul>	<p>-10 °C to +65 °C / +14 °F to 149 °F</p> <p>700 – 1060 hPa</p> <p>0 % – 95 %, no condensation</p>
Ambient conditions, transport <ul style="list-style-type: none"> <li>• Temperature</li> <li>• Air pressure</li> <li>• Humidity</li> </ul>	<p>-10 °C to +65 °C / +14 °F to 149 °F</p> <p>700 – 1060 hPa</p> <p>0 % – 95 %, no condensation</p>

Power supply (variant-dependent)	
<ul style="list-style-type: none"> <li>• Batteries <ul style="list-style-type: none"> <li>– Supply voltage</li> <li>– Battery type</li> </ul> </li> <li>• Power supply unit <ul style="list-style-type: none"> <li>– Supply voltage</li> <li>– Maximum current consumption</li> <li>– Insulated device, protection class II (EN 60601-1)</li> </ul> </li> </ul>	9 V 6 x type AA/1.5 V  12 V Typically 500 mA  <input type="checkbox"/>
Mains voltage	100 V – 240 V
Power supply frequency	50 Hz – 60 Hz
Current consumption	Approx. 55 mA
<ul style="list-style-type: none"> <li>• <b>seca 787</b></li> <li>• <b>seca 797</b> <ul style="list-style-type: none"> <li>– WiFi function and USB interface deactivated</li> <li>– WiFi function and USB interface activated</li> </ul> </li> </ul>	Approx. 55 mA Approx. 55 mA Approx. 250 mA
Maximum runtime in battery mode	Approx. 32 hours
<ul style="list-style-type: none"> <li>• <b>seca 787</b></li> <li>• <b>seca 797</b> <ul style="list-style-type: none"> <li>– WiFi function and USB interface deactivated</li> <li>– WiFi function and USB interface activated</li> </ul> </li> </ul>	Approx. 32 hours Approx. 32 hours < 5 hours, power supply unit recommended
USB interface ( <b>seca 797</b> )	
<ul style="list-style-type: none"> <li>• Current consumption of barcode scanner (optional accessory)</li> </ul>	max. 500 mA
Metrology data, measuring rod	
<ul style="list-style-type: none"> <li>• Measuring range 1</li> <li>• Measuring range 2</li> <li>• Graduation</li> <li>• Accuracy</li> </ul>	10 cm – 122 cm (4 inch – 48 inch) 122 cm – 230 cm (48 inch – 90.5 inch) 1 mm (1/8 inch) ± 5 mm
Medical device in accordance with Regulation (EU) 2017/745	Class I with measuring function
EN 60601-1: Medical electrical equipment, type B	
Type of protection	IP20
Duty cycle	Continuous duty
Application parts in line with IEC 60601-1	Weighing platform, measuring flap

## 11.2 Weight measurement

Maximum capacity	250 kg/550 lbs
Minimum capacity	0.2 kg/0.4 lbs
Graduation	0.1 kg/0.2 lbs
Tare range	up to 250 kg (subtractive)
Accuracy	± 100 g (0.22 lbs)
<ul style="list-style-type: none"> <li>• 0 kg to 70 kg (0 lbs to 132 lbs)</li> <li>• 70 kg to 250 kg (132 lbs to 551 lbs)</li> </ul>	± 0.15 %



## 12 OPTIONAL ACCESSORIES AND SPARE PARTS (SECA 787)

Accessory/spare part	Article number
Switch-mode power supply unit: 100–240 V~ / 50–60 Hz, 12 V=/0.5 A	68 32 10 270 009
<b>seca 477</b> handrail (→ <a href="#">Compatible seca products, page 47</a> )	477 0000 009
Integration (→ <a href="#">Compatible seca products, page 47</a> ) <ul style="list-style-type: none"> <li>• <b>seca 452</b> external interface module</li> <li>• <b>seca connect 103</b> software</li> <li>• Barcode scanner (medical device)</li> </ul>	452 0060 009 Will be provided in the course of integration projects Cannot be ordered from seca, for recommendations see system instructions for use of <b>seca 103/452</b>

## 13 OPTIONAL ACCESSORIES AND SPARE PARTS (SECA 797)

Accessory/spare part	Article number
Switch-mode power supply unit: 100–240 V~ / 50–60 Hz, 12 V=/0.5 A	68 32 10 270 009
<b>seca 477</b> handrail (→ <a href="#">Compatible seca products, page 47</a> )	477 0000 009
Integration (→ <a href="#">Compatible seca products, page 47</a> ) <ul style="list-style-type: none"> <li>• <b>seca connect 103</b> software</li> <li>• Barcode scanner (medical device)</li> </ul>	Will be provided in the course of integration projects Cannot be ordered from seca, for recommendations see system instructions for use of <b>seca 103/452</b>

## 14 COMPATIBLE SECA PRODUCTS

Measuring station	Handrail	Integration	
		External interface module	Software
			
787 1721 009 787 1821 009	–	<b>seca 452</b>	<b>seca connect 103</b>
787 1721 004 787 1821 004	<b>seca 477</b>	<b>seca 452</b>	<b>seca connect 103</b>
797 1721 009 797 1821 009	–	–	<b>seca connect 103</b>
797 1721 004 797 1821 004	<b>seca 477</b>	–	<b>seca connect 103</b>

## 15 DISPOSAL

### 15.1 Disposing of the device

---

Do not dispose of the device in your household waste. The device must be properly disposed of as electronic scrap. Follow your respective national regulations. For more information, please contact seca Service or an authorized service partner.

### 15.2 Disposing of batteries and rechargeable batteries

---

Do not dispose of used batteries and rechargeable batteries in household waste, regardless of whether they contain harmful substances or not. As a consumer, you are legally obliged to dispose of batteries and rechargeable batteries via local authority collection points or trade collection points. Only dispose of batteries and rechargeable batteries once they are completely discharged.

## 16 WARRANTY

There is a two-year warranty period from delivery for defects attributable to poor materials or workmanship. All movable parts, e.g. batteries, cables, power supply units, rechargeable batteries etc. are exempt. Defects which come under the warranty will be repaired for the customer free of charge against proof of purchase. Additional claims cannot be considered. Costs of transport to and from seca are the responsibility of the customer if the device is located somewhere other than the customer's headquarters. In the event of transport damage, claims under warranty can only be made if the complete original packaging was used for transport and the device was secured and fastened in it according to its originally packaged condition. You should therefore keep all packaging parts.

The warranty will be voided if the device is opened by persons not expressly authorized by seca to do so.

In the event of a warranty issue, please contact your local seca office or the dealer from whom you ordered the product.

## 17 DECLARATIONS OF CONFORMITY

### 17.1 Europe

---

seca gmbh & co. kg hereby declares that the product complies with the terms of the applicable European directives and regulations. The unabridged declaration of conformity can be found at [www.seca.com](http://www.seca.com).

## 17.2 USA and Canada

seca 797

WiFi Module:  
FCC ID: 2AC7Z-ESPWROOM02  
IC ID: 21098-ESPWROOM02

This device complies with Part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions. (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

### NOTE

This device complies with Part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions:

- This device may not cause harmful interference.
- This device must accept any interference received, including interference that may cause undesired operation.

### NOTE

Changes or modifications made to this equipment not expressly approved by seca may void the FCC authorization to operate this equipment.

### NOTE

Radiofrequency radiation exposure information: This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance of 1 m between the radiator and your body. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

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