

seca 676

Digital wheelchair scale with handrail and wireless transmission

EMR integrated



- Generously sized platform.
- Railing provides a stable standing aid.
- Foldable and mobile.
- Two fold-away ramps included.
- Pre-TARE function with three memory cells.
- At the press of a button measured results can be transmitted wirelessly to a seca 360° wireless digital printer or a PC.

seca 676:

Safe weighing – even when standing is difficult.

It is ideal to be physically fit and mobile, however, this is not possible for everyone. Many people find it increasingly more difficult to walk and stand as they get older. Many suffer from metabolic illnesses such as diabetes and often suffer from being overweight as a result. The elderly and the young suffer from injuries such as breaking a bone. Women are strongly impeded in the last months of a pregnancy. Hospitals, rehabilitation centers, and nursing homes must be prepared for all handicaps. And they will be – with the wheelchair scale seca 676.



Generously sized platform and stable railing.

Stable, functional and mobile: thanks to its large platform and sturdy handrail, the seca 676 is extremely versatile. Patients can be weighed while sitting in a wheelchair or a chair. The handrail provides valuable support for people who are frail or have difficulty walking. Its high load capacity makes the seca 676 ideal for weighing very heavy patients.

Foldable, mobile and space-saving.

After use, the scale can be folded together to save space. The sturdy locking device between the rail and platform ensures that the seca 676 stands safely even when folded. The handrail also serves as handle when the scale is folded, allowing the scale to be moved around or stowed away effortlessly and quickly on its transport castors.



Integrated ramp folds up for trouble free transport.



Intelligent functions for comfortable weighing.

The swivel-mounted display at hip level ensures comfortable operation for physicians. After being weighed, the patient can first be attended to without any stress as the weight is stored by the HOLD function. The additional weight of all aids used such as a chair or wheelchair can be deducted at the touch of a button, thanks to the TARE function. The scale has memory cells for storing three additional weights as well as the height of three persons. The BMI function provides the individual nutritional condition of the weighed person.



The scale can be folded for easy mobility.

Paperless and error-free wireless transmission of measurements.

With the seca 360° wireless technology the seca 676 can transmit measurements to the optional seca 360° wireless digital printer. With the network-capable software solutions seca analytics 115 or seca emr flash 101 and wireless USB adapter seca 456, your PC can receive and analyze measurements and forward them to an Electronic Medical Record (EMR) system. Thus, the seca 676 is EMR integrated and ready to handle digital patient records while meeting all the demands the future will bring. For more information, visit: www.seca.com

seca 676

Technical data

- Capacity: 800 lbs / 360 kg
- Graduation: 0.1 lbs / 50 g
- Dimensions (WxHxD): 36.2 x 44.1 x 45.3" / 920 x 1,120 x 1,150 mm
- Dimensions/platform (WxHxD): 31.5 x 2.2 x 38" / 800 x 55 x 965 mm
- Weight: 83.8 lbs / 38 kg
- Power supply: Power adapter
- Data transmission: seca 360° wireless technology
- Functions: TARE, pre-TARE, HOLD, auto-HOLD, BMI, auto-CLEAR, lbs/kg switch-over, damping, SEND/PRINT
- System compatible with: 360° wireless digital printer advanced seca 466, 360° wireless digital printer seca 465, 360° wireless USB adapter seca 456, software seca analytics 115 and seca emr flash 101, stadiometer seca 274 or seca 264
- RS232 compatible