

# Digital flat scale with high-quality two-component rubber surface



# seca 803

- + Two-component rubber coating offers higher security and comfort while standing on the scale
- + Four stable rubber feet provide a stable feeling
- + Large, easy-to-read 1.1 inch LCD digits
- + Optionally available in black or white
- + More than 180 years of experience for maximum precision
- + From the global market leader for medical weighing and measuring





## Start weighing medically

Go ahead, pick up our starter scale and take a close look. You'll be impressed by its classic design and its surface gleaming with high-quality materials. The high-quality rubber coating consists of a special combination of two components which lends softness and solidity. The material can stand up to all commercially available cleaning and disinfecting products

| Capacity           | 330 lbs / 150 kg  |
|--------------------|---|
| Graduation         | 0.2 lb / 100 g  |
| Dimensions (WxHxD) | 12.4 x 1.5 x 12.8 in / 316 x 37 x 326 mm  |
| Platform (WxHxD)   | 11.6 x 1.4 x 12 in / 295 x 35 x 305 mm  |
| Net weight         | 3.3 lbs / 1.5 kg  |
| Power supply       | Batteries   |
| Functions          | Auto-HOLD, Tap-on automatic switch-on, Automatic switch-off, lbs/kg switch-over |



### Sturdy design for precise results

Whether black or white, for professional or private use, this seca starter scale satisfies the highest demands for private use and anyone who appreciates high-quality workmanship. The seca 803 has its strengths on the outside and on the inside. Underneath the weighing platform you'll see longitudinal and cross braces set one centimeter apart, giving the scale unparalleled stability.



### Lightening the load

Thanks to its low tare weight of only 3 lbs and the handle built on the back, the scale is easy to carry and makes it suitable for users on the move. Viewing the measurements is just as easy from the unconventionally large display with 1.1 inch-high figure. No corners are cut with the seca 803.

We reserve the right to make modifications.