

Flat scale with foot switches and double display


seca 874

- + Double display lets the doctor and patient read simultaneously
- + Reliable and durable: seca quality that pays off
- + May be used standing up thanks to the foot switches
- + Can weigh 10,000 times with a single set of batteries
- + Mother+child function to weigh toddlers



seca 874

Flat scale with foot switches and double display



Setting up, tapping, weighing

Equipped with leveling feet and a level, the seca 874 flat scale can be safely set up horizontally even on uneven surfaces. One battery pack is sufficient for thousands of weighings and makes the seca 874, together with a carrying handle and low weight, well-suited for mobile use. Tapping the foot pedal activates the scale, the other foot pedal operates the mother+child function, allowing a child to be weighed on the mother's arm.

Read twice as well

A double display has the decisive advantage that it can be read from both sides. The doctor can therefore quickly record the result, regardless of whether he is in front of or behind the patient. And the patient feels involved, because the scale also shows them their weight.



Capacity	440 lbs / 200 kg
Graduation	0.1 lb < 330 lbs > 0.2 lb / 50 g < 150 kg > 100 g
Dimensions (WxHxD)	12.6 x 2.4 x 14.2 in / 321 x 60 x 361 mm
Platform (WxHxD)	11.3 x 2.4 x 11 in / 288 x 60 x 280 mm
Net weight	9 lbs / 4.1 kg
Power supply	Batteries
Functions	Auto-HOLD, Automatic switch-off, lbs/kg switch-over, Mother+child function
Accessories	+ Stadiometer seca 217 in combination with adapter element seca 437 + Carrying cases seca 409 / seca 415 or seca 414 (for combination with stadiometer and adapter element)

Quality from seca: sturdy and rainforest tested

In elaborate testing, not only is durability tested, but also resistance to humidity and heat as well as impact resistance. A stable steel frame, ribbed components and high-quality ABS plastic, from which motorcycle helmets are also manufactured, ensure that seca scales meet all requirements. A long scale life is thus guaranteed, whether in an air-conditioned room or in a temporary practice in the rainforest.

We reserve the right to make modifications.

seca north america | west
13601 Benson Avenue
Chino, CA 91710
USA

seca north america | east
1439 Stuart Engals Blvd, Suite 200
Mount Pleasant, SC 29464
USA

T +1 800 542 7322
F +1 888 705 7397
E info.us@seca.com

Precision for health
seca.com