

Flat scale with Mother+child function

seca
Precision for health


seca 876

- + Precise weighing technology for medical use
- + Mother+child function to weigh toddlers while cradled on the mother's or father's arm
- + May be switched on while standing up thanks to the tap-on function
- + Leveling feet for strong stability
- + Can weigh 10,000 times with a single set of batteries



seca 876

Flat scale with Mother+child function



The classic flat scale for medical use

The seca 876 is designed to deliver the most precise results. With its steel frame, it is extremely stable and at the same time light enough to be used on the move. The knotted rubber surface makes the scale slip-proof and particularly safe.



Setting up, tapping, weighing

Equipped with levelling feet and a level the seca 876 flat scale can be safely set up horizontally even on uneven surfaces. One battery set is sufficient for thousands of weighings and makes the seca 876 independent of the mains power supply. Tapping the tray activates the scale. A touch of the button activates the Mother+child function to allow a child being weighed while cradled on the mother's arm.

Capacity	550 lbs / 250 kg
Graduation	0.2 lb < 330 lbs > 0.5 lb / 100 g < 150 kg > 200 g
Dimensions (WxHxD)	12.6 x 2.4 x 14 in / 321 x 60 x 356 mm
Platform (WxHxD)	11.3 x 2.4 x 10.9 in / 288 x 60 x 278 mm
Net weight	8.6 lbs / 3.9 kg
Power supply	Batteries
Functions	Tap-on automatic switch-on, Automatic switch-off, lbs/kg switch-over, Mother+child function
Accessories	+ Stadiometer seca 217 in combination with adapter element seca 437 + Carrying cases seca 409 / seca 415 or seca 414 (for combination with stadiometer and adapter element)



Quality from seca: sturdy and rainforest tested

In elaborate testing, not only is durability tested, but also resistance to humidity and heat as well as impact resistance. A stable steel frame, ribbed components and high-quality ABS plastic, from which motorcycle helmets are also manufactured, ensure that seca scales meet all requirements. A long scale life is thus guaranteed, whether in an air-conditioned room or in a temporary practice in the rainforest.

We reserve the right to make modifications.

seca north america | west
13601 Benson Avenue
Chino, CA 91710
USA

seca north america | east
1439 Stuart Engals Blvd, Suite 200
Mount Pleasant, SC 29464
USA

T +1 800 542 7322
F +1 888 705 7397
E info.us@seca.com

Precision for health
seca.com

5/2023